



# dbLabCal

## Manual



**dbLabCal**  
The LIMS for bioanalytics

- Analytical results management
  - Data import from HPLC,GC+LC/MS-MS,Immunoassays
- Calculations and statistical evaluations
- Acceptance checks
- Chromatographic data evaluation
- QC, QA, ES according CRF21 Part11, Audit Trail

*Copyright © 1994-2015 Milan Vagaday*

**Index of Contents:**

<b>INTRODUCTION .....</b>	<b>4</b>
PROGRAM TASKS .....	6
TYPICAL WORKFLOW .....	7
USER AUTHORIZATION .....	8
NAMING CONVENTIONS .....	9
<b>OPERATION .....</b>	<b>11</b>
GENERAL .....	11
<i>Keyboard and Mouse</i> .....	11
<i>Screen</i> .....	13
LOGIN .....	14
PROJECT MENU .....	15
<i>New</i> .....	15
<i>Load</i> .....	18
<i>Edit</i> .....	20
<i>Options Menu</i> .....	21
Print/Export Options .....	21
Graphics .....	23
Number Format .....	23
SUB Format .....	24
Display Time as .....	24
Display LOQ as .....	24
Display Not Correct SUB results as .....	24
Acceptance Criteria .....	25
Reporting Texts .....	27
Set Target Folders .....	28
Standard Column Width! .....	28
<i>Analyte</i> .....	29
<i>Regression Model, Weighting, Readings</i> .....	29
<i>Audit Trail and Electronic Signature</i> .....	30
<i>Print</i> .....	31
<i>Export</i> .....	33
<i>Changing Project Status</i> .....	35
<i>Results for HoLaRo... / Results for BI</i> .....	36
<i>Results in ASCII File!</i> .....	36
<i>Exit</i> .....	36
RESULTS MENU .....	37
<i>Statistics</i> .....	38
<i>Subject Samples Results</i> .....	40
<i>Validation Samples</i> .....	44
<i>Re-Assays</i> .....	50
Samples to be Re-analyzed .....	50
Re-assayed Samples .....	51
Incurred Samples .....	55
<i>Excluded Values</i> .....	56
<i>Date of Analyses</i> .....	56
BATCH MENU .....	56
<i>Import File</i> .....	57
<i>Import from Empower2</i> .....	60
<i>Load</i> .....	61
<i>Print</i> .....	62
<i>Export</i> .....	63
<i>Change Batch Status</i> .....	64

---

BATCH DATA MENU.....	65
<i>Batch List</i> .....	67
Editing Sample Names.....	68
Editing Chromatogram Flags.....	69
CHROMATOGRAPHY MENU .....	71
DB MENU.....	76
MENÜ EXTRA.....	83
<b>TIPS AND TRICKS.....</b>	<b>84</b>
<b>ERRORS AND ERROR MESSAGES.....</b>	<b>90</b>
<b>APPENDIX 1: HOLARO- AND. BI- ASCII-FILES.....</b>	<b>91</b>
<b>APPENDIX 2: PE SCIEX LC-MS (ANALYST).....</b>	<b>94</b>
<b>APPENDIX 3: EMPOWER2.....</b>	<b>95</b>
<b>APPENDIX 4: XCALIBUR.....</b>	<b>97</b>
<b>APPENDIX 5: SOFTMAX PRO.....</b>	<b>98</b>
<b>APPENDIX 6: MAGELLAN (ELISA READER).....</b>	<b>100</b>
<b>APPENDIX 7: ACCESS2.....</b>	<b>101</b>
<b>APPENDIX 8: ISE.....</b>	<b>102</b>
<b>APPENDIX 9: ICP-MS ELAN.....</b>	<b>103</b>
<b>APPENDIX 10: ICP-HPLC-MS CHROMERA.....</b>	<b>108</b>
<b>APPENDIX 11: FACS.....</b>	<b>109</b>
<b>APPENDIX 12: SEARCHLIGHT.....</b>	<b>109</b>
<b>APPENDIX 13: MESOSCALE (MSD).....</b>	<b>110</b>
<b>APPENDIX 14: CALCULATION OF WEIGHTED LINEAR REGRESSION.....</b>	<b>111</b>
<b>APPENDIX 15: FLOWCHART: REPORTING FINAL RESULTS (LANG&amp;BOLTON).....</b>	<b>112</b>
<b>APPENDIX 16: FLOWCHART: REPORTING FINAL RESULTS (HOLARO).....</b>	<b>113</b>

## **Introduction**

The aim of DBLABCAL is the management of analytical data and results generated during bioanalytical studies (e.g. bioavailability, bioequivalence study). It is compatible with a wide variety of analytical methods and instruments (HPLC, LC/MS, GC/MS etc.).

DBLABCAL is an interface to a database in which the information on analytical studies, on the batches of the studies and, lastly, "all" data of the individual chromatograms are stored. DBLABCAL sends the user's commands or questions (selected from the individual menu items) to the database and then displays/prints/exports the answers from the database in a formatted and meaningfully arranged way. In addition, many plausibility tests are performed to support the user.

All data of a chromatographic run (batch) are imported into the database by an ASCII file.

DBLABCAL monitors the data for any irregularities and provides information/hints about certain situations. Furthermore, the program contains user management functions and only actions which are meaningful at that particular time and for which the user has access authority can be performed

All actions (commands) of the user which modify the data/results of a project are documented in the audit trail. DBLABCAL fulfils the CFR21Part11 requirements.

The chromatographic source data (retention times, peak heights, peak areas, calculated concentration, etc.) cannot be modified with DBLABCAL once they have been imported.

It is also possible in DBLABCAL, based on the analytical data stored in the DBLABCAL database, to extract management information data such as project status or number of samples analyzed in a specific period of time in a specific department, etc.

In order to work with the database the user must comply with certain rules for assigning sample names. Then, all data of the project are available immediately, free of errors, for further processing (e.g. biometric evaluation, export into other programs for report generation or for further statistical analyses, etc.).

This version of DBLABCAL works with all Oracle database versions. Versions from 8i to 11gR2 as well as Oracle XE versions were tested.

DBLABCAL works in two different languages. The user can freely switch between the two. Furthermore, the administrator can modify and adapt all annotations, messages, texts and system texts at any time.

This manual describes only the most important points. There may therefore be some situations in which the program returns questions or messages to the user that are not explicitly described herein.

## **Program Tasks**

The major tasks of DBLABCAL can be summarized as follows:

### **Administration:**

- Project Administration
- Permission Administration
- Display, printing and export of project results
- Export of subject results into ASCII files for further processing
- Release of project results

- Import of chromatographic batch data
- Display, printing and export of data from individual batches
- Release of individual batches

- Editing of sample names
- Release of individual chromatograms

- Long Term Stability Planning

### **Calculations:**

- Calibration curves
- Concentration of unknown samples  
(QCs, subject samples, validation samples, etc. ...)
- Statistical evaluation of project results

### **Monitoring of all results for correctness:**

- Verifying the correctness of project results
- Verification that batch data meet acceptance criteria

### **Chromatography data analysis:**

- Display of retention times and peak widths for peaks of all project chromatograms  
e.g. for verifying the correctness of integration, peak assignment, etc.
- Selection of specific data from the database  
e.g. peak areas of all QCs to demonstrate their stability throughout the duration of measurements...

### **Audit trail / Electronic Signature**

### **Management information**

## Typical Workflow

DBLABCAL is generally used as follows:

	Person	Action	Comment
<b>BEFORE SAMPLE MEASUREMENT</b>			
1	PI/Study Director (Analyst)	defines new project	Define project code, comment, access authority for analyst, analyte, IS name(s), matrix, parameters, model, weighting, conc. unit(s) analytical method(s), project type, if routine project: number of samples to be analyzed  see page 15
<b>DURING SAMPLE MEASUREMENT</b>			
2	Analyst	analyzes batches	Writes batch lists Runs batches
3	Analyst	checks chromatograms BEFORE data export into dbLabCal	Checks chromatography data for analyte(s) and ISs Sets the acceptance on the chromatograms (flags)
4	Analyst	creates export ASCII files (*.lca, *.xls, *.asc, *.csv...)	Using the chromatography software (Analyst, Xcalibur, Empower Magellan, etc.)
5	Analyst	imports ASCII files into DBLABCAL	  see page 57
6	Analyst	sets chromatogram flags in dbLabCal if necessary  <i>Chromatogram release</i>	In accordance with the notes on the chromatograms  see page 67
7	Analyst	sets batch status (to "batch accepted" or "not accepted")  <i>Batch release</i>	depending on the acceptance criteria check results  see page 64
8	PI/Study Director	reviews batch data<->chromatograms	informs the analyst if necessary (possible also via the batch memo) sets R / C flag if required
continue at 2			
<b>AFTER SAMPLE MEASUREMENT</b>			
9	Analyst	sets project status to 'finished'  <i>Project release</i>	  see page 34
10	PI/Study Director	sets project status to 'released' (it also locks all batches, batches remain locked even after project status is reset) <i>Project release/lock</i>	  see page 34
11	QC-Check	QC check of the whole project	Consistency and correctness re. project plan, chromatographic data etc...
12	QM/QA	Documents the quality audit in the QM/QA department	
13	Any user	prints, exports batch data prints, exports project results	as required for documentation and report

Data are released/locked in three steps:

- release of individual chromatograms (flag)
- release of individual batches (batch status), "batch lock"
- release of whole project (project status), "project lock"

## User Authorization

The following permission levels exist:

ReadOnly	No permission to edit data
Analyst	may modify data of "own" projects
PI/Study Director	PI/Study director may modify data of "own" projects
QC/QA Reviewer	permission to document QC/QA work otherwise like ReadOnly
Department Manager	PI/Study director permissions for all projects in the department
Administrator	permissions for all projects

Action	ReadOnly	Analyst	PI/Study Director	QC/QA	Comment
Create new project		+	+		Analyst cannot assign other user
Edit project data		+	+		if project "released" add new analyte(s) only allowed
Grant permissions for project			+		
Change model, weighting and readings		+	+		if project type "validation"
Change acceptance criteria			+		as long as project not "released"
Change settings and output texts	(+)	+	+		ReadOnly-User only for current session, save not possible
Print/export of data	+	+	+	+	
Finish project		+	+		
Release project, Reset project status			+		
Data export in ASCII			+		if project "released"
Import batch		+	+		as long as project not "released"
Change unit, batch number, start date		+	+		Analyst only if project „started“ PI/study director only as long not "released"
Change batch status		+	+		Analyst only if project „started“ AND to a batch status in accordance with batch acceptance criteria
Final result selection (from re-assays)		+	+		Analyst only as long as project „started“
Change chromatogram flag		+	+		Analyst only as long as batch status "not set yet"
Change sample name		+	+		Analyst only as long as batch status "not set yet"
QC check, QA				+	after project was "released"

+ permitted



## Naming Conventions

### Batches:

A batch file (ASCII file for import into DBLABCAL) can have any name. The file extension describes the ASCII file format

Extension	*.lca	created by Analyst (Sciex)
	*.xls	created by Xcalibur, LabX, MultiLab Pilot, SearchLight
	*.rep	created by Elan-ICP
	*.csv	created by Access2, FACS, Mesoscale (MSD)
	*.asc	created by Magellan (ELISA)
	*.txt	created by SoftMax Pro

It is recommended to use ASCII file names consisting of any 2 letters and 2 numbers which define the batch number, e.g. aa01.lca, mn17.xls, xy99.rep, de5an01.csv and de5an02.asc etc.

### Samples:

Sample Type	Code	Further Information	Example																	
System Suitability Sample	<b>SSS</b>		<b>SSS</b>																	
Standards	<b>CAL</b>	Nominal concentration	<b>CAL</b> 10.0 ng/mL <b>CAL</b> 10.0/25.0 ng/mL																	
QC's	<b>QCS</b>	Nominal concentration	<b>QCS</b> 150 ng/mL <b>QCS</b> 150/450 ng/mL																	
Unknown samples	<b>SUB</b>	Number, Period, Sampling time (ALL VALUES MUST BE NUMERIC!)	<b>SUB01DG01TP2d4h30m</b> <b>SUB15DG02TP2.5h</b> <b>SUB99DG03TP5</b>																	
		Number or name	<b>SUB123456</b> <b>SUBMD44_F01</b>																	
Validation samples	<b>VAL</b>	Nominal concentration																		
		<table border="1"> <thead> <tr> <th>"Matrix"-Group *)</th> <th>"Temp"-Group *)</th> <th>Duration</th> </tr> </thead> <tbody> <tr> <td><b>E</b> extracts</td> <td><b>R</b> room temp.</td> <td>given in hours or _d_h_m format</td> </tr> <tr> <td><b>N</b> matrix</td> <td><b>K</b> refrigerator</td> <td></td> </tr> <tr> <td><b>P</b> thawing cycles</td> <td><b>G</b> freezer</td> <td></td> </tr> <tr> <td><b>A</b> other (whole blood)</td> <td><b>T</b> deep-freezer</td> <td></td> </tr> <tr> <td><b>B</b> n/a</td> <td></td> <td></td> </tr> </tbody> </table>	"Matrix"-Group *)	"Temp"-Group *)	Duration	<b>E</b> extracts	<b>R</b> room temp.	given in hours or _d_h_m format	<b>N</b> matrix	<b>K</b> refrigerator		<b>P</b> thawing cycles	<b>G</b> freezer		<b>A</b> other (whole blood)	<b>T</b> deep-freezer		<b>B</b> n/a		
	"Matrix"-Group *)	"Temp"-Group *)	Duration																	
<b>E</b> extracts	<b>R</b> room temp.	given in hours or _d_h_m format																		
<b>N</b> matrix	<b>K</b> refrigerator																			
<b>P</b> thawing cycles	<b>G</b> freezer																			
<b>A</b> other (whole blood)	<b>T</b> deep-freezer																			
<b>B</b> n/a																				
<b>QCS</b>	Nominal concentration and dilution factor (df)	<b>QCS</b> 1000 ng/mL df in the chromatography software																		

\*) The codes for "Matrix" group (E, N, P, A) or "Temp" group (R, K, G, T), respectively, are assigned by default to the matrix type and storage conditions as shown in table.

It is possible to assign any description to any "Matrix"- "Temp" combination.

E.g., VAL **AT** 2h xxx pg/mL may belong to the validation/stability samples group:  
„Plasma stabilized with X and kept cool in ice water during sample collection”.

It is required for CAL, QCS, and VAL samples to enter the nominal concentration after one or more blanks. Nominal concentrations should be entered in the order of their retention times separated by '/' in case the CAL, QCS, or VAL sample contains more than one analyte with different nominal concentrations.

Examples:

CAL 500 ng/mL, VAL NG48 10/20/20 ng/mL, QCS 5/10 ng/mL etc. If all analytes are of the same nominal concentration it is sufficient to enter the nominal concentration just once.

**dbLabCal read nominal concentrations from Analyst, Empower2 and Xcalibur DIRECTLY. Nominal concentrations in the sample name are ignored for Analyst, Empower2 and Xcalibur. They are used ONLY if this information was missing in the chromatography software.**

In the case of SUB samples, it must be decided before starting the measurements whether the samples should be named SUBxxxxDGxxTPxxx or only SUBxxxxx as both modes are not allowed simultaneously in one project.

It is possible to switch between the SUB naming modes (SUBxxxxDGxxTPxxx or SUBxxxxx). Potentially missing Period and Sampling Time data may need to be completed.

The SUBxxxxDGxxTPxxx naming is appropriate only when the samples are indeed samples from a bioequivalence study with "real" periods and time points.

If, for example, samples from a Phase II study or a toxicological study are to be analyzed, the sample names will probably better be coded in such a way that the results appear in RESULTS / SUBJECT SAMPLES in the desired order (the results are sorted alphabetically) and consequently, can easily be processed further (e.g. in Excel).

**Subject Samples Results are sorted in alphabetical order. That means that also the numbers are sorted in alphabetical order. (example: 1, 10, 11, 2, 222, 3 etc.).**

**In case the samples should be sorted as: 1, 2, 3, 10, 11, 222 etc, it is necessary to write numbers in SUBxxxxxx as SUB0001, SUB0002, SUB0003, SUB0010, SUB0011, SUB0222, etc.**

For instance, a sample from a female rat no. 118, dosage 10 mg/kg i.v. obtained on day 7, 12 hours, could be identified as follows: SUBIV\_D010\_F\_D07H12\_NO0118 (or short, e.g., SUBI\_010\_F\_07\_12\_118). The results will then be shown sorted by type of dosage, dosage, sex, time and animal number. If it were necessary to separately report the results first by sex and then by dosage, the samples would be identified as follows:

SUBF\_IV\_D010\_D07H12\_NO0118. The separator '\_' is useful for improving readability and facilitating the later extraction of individual data with SEARCH AND REPLACE into Excel where they will be written into separate columns.

The determining factor for naming samples is the way of subsequent result data processing and reporting.

**Operation**

**General  
Keyboard and Mouse**



Both keyboard and mouse operations follow the Windows conventions. Menus, availability of certain commands, etc. are adapted dynamically to the current data and the access authority of the user. Furthermore, certain menus and commands are disabled if their use is not meaningful or not allowed at that particular point in time. Whenever the OK and CANCEL buttons are shown on the screen, the ENTER and ESCAPE keys are available for use. Pressing the ENTER key is equivalent to clicking on the OK button and, similarly, pressing the ESCAPE key has the same effect as clicking on the CANCEL button.

To terminate an input and jump to the next element (e.g. to the next text field of the current dialog), the TAB key and not the ENTER key is used in Windows.

The mouse is used in very much the same way as it is in Windows:

- Left mouse button: Click: Select
- Double-click: Select with confirmation (OK)
- Right mouse button: Click: Selection menu / context menu, where available

Some inputs can be performed by more than one way. For example, the project comment can be entered via the menu PROJECT | NEW... or PROJECT | EDIT... or by double-clicking on COMMENT ON THE PROJECT.



dbLabCal accept following shortcuts:

ENTER/ESCAPE	Default keys for clicking on OK or CANCEL buttons, respectively
CTRL-INS or CTRL-C SHFT-INS or CTRL-V	In addition to the menu-driven export (see below) it is possible at any time to export (marked) data and/or graphics from the active window into other programs via the Windows clipboard (as generally done in Windows).
STRG+L, +P, +E	Laden, Drucken, Export des Projekts, bzw. der aktuellen Sequenz
CTRL - +/- (NUMERIC KEYPAD)	loads next / previous batch
SHFT+CTRL+C, +Q, +V, +S, +K, +D	Sample type filter on/off (in batch view)
SHFT+CTRL +R	Current batch is immediately re-calculated
CTRL-R	(un)marks „Random Repeats“/IncurredSamples ISR are marked by yellow background color
press SHIFT while changing the concentration unit in PROJECT   EDIT...	nominal concentrations of all until now analyzed CAL, QC and VAL samples will be re-calculated according new concentration unit. (factor 1000, 1 000 000 etc.)
press CTRL while changing the concentration unit, readings or regression model in PROJECT   EDIT...	changes parameter for CURRENT ANALYTE ONLY (default behavior is to do changes the parameter for all analytes in one chromatographic run)
press CTRL at program start...	Project load dialog will always be displayed
press CTRL while clicking on the user name (status line, right side)...	Login dialog will be displayed for re-login
press CTRL while changing chromatogram flag...	chromatogram flag will instantly be changed for all analytes of the chromatogram
press CTRL while double clicking on a value or batch number in any results table or in the chromatography menue... (or press CTRL+G)	dbLabCal “jumps” to the chosen batch and chromatogram return to original view by pressing the ESCAPE key
CTRL while certain actions gives the PI/study director more freedom than normally available, e.g. before selecting final results from re-assays (the so-called ‘study director special wish key’)	

## Screen

Title bar: Shows the current information for  
program name / version - project code / analyte name / batch number

Menus:

Status line: The status line shows the following information:

left: progress display or time

center: - messages to current actions of the program

- explanations to the significance of the fields under the mouse pointer,

and

- tips to possible user actions

- name of current data base

right: name of current user


Between the menus and the status line the screen is hierarchically divided into three segments:


top: information to current project and to analyte

center: information to current batch (if batch data view was chosen)

bottom: data of current batch in the selected view  
or data (results) of current project

The MEMO fields allow for entering comments or notes to results or to the individual batches. They are also suitable vehicles information exchange. The MEMO field displays green if there is content, otherwise they are dark gray.

Occasionally the button  is displayed. Clicking on this button opens a dialog offering more settings available for the current display.




	ng/mL	p.h.t.	ng/mL[calc]	dev [%]
1	2.50	0.0340	2.47	-1.1
2	6.25	0.0859	6.41	2.5
3	12.5	0.1645	12.4	-1.1
4	25.0	0.3302	24.9	-0.2
5	62.5	0.8203	62.1	-0.6
6	125	1.6509	125	0.1
7	250	3.2914	250	-0.2
8	625	8.3357	632	1.2
9	1250	16.3935	1240	-0.5
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				

---

## Login

At program start the following Login Dialog appears:




The image shows a standard Windows-style dialog box. The title bar is blue and contains the text "LOGIN (dblabcal-db/dblabcaldev)". The main area of the dialog is white and contains two labels, "UserID" and "Password", each followed by a text input field. At the bottom of the dialog, there are two buttons: "EXIT" on the left and "LOGIN" on the right.

User can log into DBLABCAL even if another user is logged into Windows. The user is locked after a specific inactivity time period.


## Project Menu

Studies are managed through this menu



## New...

A new dialog for entering all required information for a new project appears:



- Project code

*By default a project code is allowed just once (per department). In certain situations, however, it may be advantageous to subdivide a project into, for instance, a validation part and a routine part, or subdivide by matrix etc. In such a case the department manager, after acknowledging the security prompt, can create another project with an existing project code in the data base.*

- GLP code for the project

- Sponsor code for the project


- Comment to project - usually the analyte(s) name and the matrix

- Department, analyst and PI/study director, deputy of the PI

- Analyte names measured in the current project, plus other details such as project types, analytical method, matrix, readings, concentration unit, model, weighting, and, if used, name of internal standard and decision whether the program will or will not set the '+' flag in chromatograms of pre-dose samples. (see page 67 for more information about flags)

If the current project is a "routine" project the number of samples to be measured must be entered in the last column and the used SUB naming type can be chosen

User can select in last column (right mouse) if long term stability is due for the selected analyte/matrix. He/she will be reminded 2 weeks before the due date to perform the stability evaluations.

It is possible to copy all analytes and their settings from another project by clicking  on next to PROJECT-No. into current project.

Furtermore, it is possible with COPY ALL SETTINGS FROM: to copy all settings (submenu PROJECT | SETTINGS) from another project into the current project.

If several analytes are measured in a project not all analytes in a chromatogram have to be determined. It is thus possible, for instance, to measure analyte X and metabolites 1, 2 and 3 in plasma in one 'run' with analyte X in urine and, perhaps, analyte X (total) in urine being part of the project XY001 as well. This example would lead to 6 analyte entries.

*If a project includes several chromatographic methods and/or different matrices it is more reasonable to create a separate projects (PI or administrator are allowed to create dbLabCal projects with identical project codes) for each of the chromatographic method and/or matrix in DBLABCAL.*

A list with all analytes stored in dbLabCal appears after double-click in the text field ANALYTE or INT.STD. Double-click on analyte name select that analyte for the new project. This saves typing work as well as using exact the same names for same analytes makes any further evaluation of analyte data beyond 'project limits' easier.

When entering an analyte the fields for internal standard, matrix, readings, model, weighting etc. will automatically be filled with the contents of the upper line.



By clicking on >>> in the list of analysts the PI/study director can name up to 10 analysts (use SHFT or CTRL key) who will also be granted access to the newly created project. By clicking on the <<< list of analysts will be closed again.

Furthermore, up to 10 user with the QC permission (ReadOnly) may be defined. Analysts and QC check user must not be identical.

If admin has set the database external access option, the user may grant or revoke access to the project also to one external user.

Department: [ ] Ext. access granted to: [NONE]  
 Study Director: [ ]  
 Analyst: [ ] >>>

An external user has to be defined first in dbLabCal through the db users dialog (page 81) in the same way as the “internal” user. Only the checkbox “external user” has to be selected. An external user may get access to the released data only or to all data.

Users  
 ID: 0 [New] [Exit]  
 User ID: EXT01  
 Name: ExtUser01  
 Password for Oracle Login: \*\*\*\*\*  
 Department: EXT  
 Authorization: Read Only  active  
 'w' flag allowed  
 Reviewer  
 external user  released data only  
 SQL Direct: No Rights  
 Capacity Browser: No Rights  
 [Find!] [Cancel] [Save!]  
 [First] [Previous] [Next] [Last]

*The program always checks all user entries for plausibility and automatically corrects non-plausible entries. If, for instance, you select 'p.a.r.' for readings without having entered internal standard before, readings is automatically corrected to 'area'. The dialog cannot be closed by selecting OK before all required data are entered.*

**Load...**

A list of all projects sorted by project code appears. 'Filter' allows selection of projects to be displayed in the project load dialog.

Each user only sees his/her own studies. All project's analytes are listed below the project list together with their status. After selecting the project and the analyte the user can load the data by clicking on OK. (or by double-click on the project or the analyte)

*Quick access to dialog: Double-click anywhere within the screen area 'Project'.*

The dialog LOAD PROJECT is also displayed at program start if automatic display of the last project data was not selected at the time of the last program exit (see section EXIT). The LOAD PROJECT dialog always displays if the CTRL key is pressed and held down during program start.

Keeping CTRL key pressed while opening the dialog let open PROJECT LOAD dialog in ReadOnly Modus. It mean that user see all projects and that loaded project also opened in ReadOnly modus. User cannot edit such projects independently of the actual users permissions.

Load Project						
Project	Comment	Project Type	changed	Sponsor Code	GLP-Code	
GA323	5-ASA / Ac-5-ASA in Human Plasma and Faeces	Routine	06-Aug-1998			BA
HA107	Ac-5-ASA / 5-ASA	Routine	27-Mrz-1998	98-0435-001-L1		BA
KA001	5-ASA and Ac-5-ASA re-validation in urine	Validation	15-Mrz-2004			BA
KA018	5-ASA/N-Acetyl-5-ASA in Plasma and Urine	Routine	09-Aug-2000	CRD-00-15		BA
KA448	5-ASA and Ac-5-ASA in plasma	Routine	15-Mrz-2001	CRD-PK-00-42		BA
KA448	5-ASA and Ac-5-ASA in urine	Routine	23-Apr-2001	CRD-PK-00-42		BA
LA131	5-ASA/Ac-5-ASA	Routine	10-Okt-2001			BA
MA188	5-ASA / Ac-5-ASA in human plasma	Routine	11-Jun-2002	PPL-559		BA
NA540	5-ASA and Ac-5-ASA in human plasma/urine	Routine	13-Okt-2004			BA
NA556	5-ASA and Ac-5-ASA in human plasma/urine	Routine	26-Apr-2011	IA358/CLAV-P-03		BA
QA435	5-ASA and Ac-5-ASA in human plasma	Routine	14-Mai-2007	TP0315		BA
TX021	5-ASA/Ac-5-ASA in Humanplasma	Validation	16-Jul-2009			BA

Analyte	Matrix	Status	Method	Project Type	Filter	
AC-5-ASA	plasma	data released	HPLC-FL	Routine	all	all
5-ASA	plasma	data released	HPLC-FL	Routine	last week	Routine
					last month	Validation
					all	all
					defined	QC in-progress
					started	QC-ed
					finished	QA-ed
					data released	Ac-5-ASA
					archived	
						Cancel
						Ok

After clicking on OK the data for the selected project and the analyte are loaded.

If another user already works on the selected project the corresponding user name is shown to the right of the analyte name. However, in this case the project data cannot be edited.

**Project data can only be edited by one user at a time.  
Other users always are in 'ReadOnly' mode.**

If the authorized user has made any modifications to the project data a message is displayed for all other users, logged in the project. After double-click on the message **New Data** the new data will be re-loaded.

---

**Edit...**

The same dialog as in the menu PROJECT | NEW is displayed.

This dialog allows for subsequent modifications to some project information, such as sponsor code, comment etc. The PI/study director is also authorized to change the project code.

The name(s) of the analyte(s) and the matrix can also be modified - the assignments of data acquired so far remains the same - or additional new analytes can be entered into the project. Existing analytes can be deleted from the list as long as the analyte status is 'applied'.

If readings, regression model and weighting are changed the complete project will be recalculated according to the new settings after exiting the dialog. Modifications to readings, regression model and weighting are automatically applied to all analytes measured in a chromatogram.

Rows for all analytes measured in a chromatogram with the current analyte (current analyte is displayed in the title bar and in the „project“ frame) have white back color. The other rows are grey.

The decision whether DBLABCAL should set or not set the '+' flag can also be changed. The complete project will be recalculated after closing the dialog.


The concentration unit can be changed as well. There are two possible scenarios. Either only the concentration unit (concentration unit text only) is changed or - if the SHFT key is pressed when selecting the new concentration unit - the nominal concentrations of the calibration, QC and validation samples measured so far can be recalculated according to the new concentration unit as well.

By default, concentration unit, readings and regression model etc. will be changed for all analytes measured in a chromatogram simultaneously. To change the parameter for one analytes only, keep the CTRL key pressed.

See also page 12 for SHIFT and CTRL keys.

All setting (all data in menu Project | Settings) can be copied from an existing project. See also page 16.

## Options Menu



Various settings applicable to the current project are selected through the options in this menu.

Acceptance criteria can only be changed by the PI/study director. Analysts and visitors are only allowed to view them but can modify other settings. However, modifications done by visitors are not saved and therefore only apply to the current session.

All setting (all data in menu Project | Settings) can be copied from an existing project in the Project | New... or Project | Edit... menu. See also page 16.

### **Print/Export Options...**

The dialog Print/Export Options... allows for changing the appearance (and the content) of the data to be printed or exported. The appearance of the printout can also be modified (frame.... header shaded..).

The number of columns/page can be set via MAX.COLUMNS (LIST OF DATA) for printing the subject results in the LIST OF DATA format for standardizing the appearance of all printouts. For instance, if 30 subjects have been measured, a good setting would be 10, or 15 (or maybe 6... try it out...).

Subject results in the STANDARD format are usually printed with 2 columns/page. However, if some of the comments to the results take up more space and extend into the next column in the printout it is better to print with 1 COLUMN ONLY (SUB, STANDARD) to generate a one-column printout. 1 COLUMN ONLY is the default setting.

**Print/Export Options**

Printer: PDFCreator (Ne01:)

Font: Arial Narrow 9

Study Code XY123

Line Spacing: 1.1

Columns:

Max. Columns (List of Data): 12

Max. Columns (CAL, QC, VAL): 12

1 Column only (SUB, Standard)

Border (Tables/Graphics)  full

Header Shaded  bold

use dec.tabstops in Word tables

Page Setup:

portrait  landscape

Page Width [cm]: 17.0

left: 2.0 top: 2.0

right: 2.0 bottom: 1.0

Column Widths:

Study: col	1	2	3	4
[cm]	2.1	9.2	3.0	2.7

Seq.: col	1	2	3	4	5	6
[cm]	2.1	7.3	1.9	3.0	0.8	1.9

Compact Presentation

set to default values

Cancel Apply Apply and Save


The settings for margins, font, font size and line spacing allow for adapting the appearance of the DBLABCAL printouts to the report. The user should "play" with these settings to find the optimum configuration.

Via HEADER COLUMN WIDTHS the user can modify the appearance (column widths) of the project and/or sequence headers on the printout and in the formatted export file to MS Word. The page width to be used can also be set for the formatted export into MS Word. Changing page width standardizes the individual column widths to the new page width.

### Graphics...

DBLABCAL adapts the scaling of all graphics to the current data by default. As desired, all or only some of the axis parameters can be set to fixed values.

These settings are accepted for all analytes in a chromatogram if they have not been individually set before. In other words, scaling can be but does not have to be the same for all analytes from a chromatogram.




The 'Graphic Options' dialog box is divided into three sections: Calibration Curve, QC Samples, and Subject Samples. Each section contains checkboxes for 'X max' and 'Y max' (or 'Y min - Y max' for QC Samples) with a default value of 'auto' and an empty text input field. The 'Subject Samples' section has checkboxes for 'X min - X max' and 'Y max' with 'auto' and empty text input fields. 'Cancel' and 'Ok' buttons are at the bottom.

Quick access: double-click on the graphic.

### Number Format...

For each reading either the number of decimal digits or the number of significant digits can be set. These settings apply to all analytes from a chromatogram.



The 'Number Format' dialog box shows a table with two columns: 'Decimal Digits' and 'Signif. Digits'. The 'Concentration' row is highlighted. Below the table are 'Cancel' and 'Apply' buttons.

	Decimal Digits	Signif. Digits
Concentration	<input type="text" value=""/>	<input type="text" value="3"/>
Time	<input type="text" value="2"/>	<input type="text" value=""/>
Dilution Factor	<input type="text" value="1"/>	<input type="text" value=""/>
Statistics [%]	<input type="text" value="1"/>	<input type="text" value=""/>
Retention Time	<input type="text" value="2"/>	<input type="text" value=""/>
Peak Area	<input type="text" value="0"/>	<input type="text" value=""/>
Peak Height	<input type="text" value="0"/>	<input type="text" value=""/>
Readings	<input type="text" value="4"/>	<input type="text" value=""/>
a	<input type="text" value="5"/>	<input type="text" value=""/>
b	<input type="text" value="5"/>	<input type="text" value=""/>
c	<input type="text" value="5"/>	<input type="text" value=""/>
r	<input type="text" value="5"/>	<input type="text" value=""/>
reassays diff%	<input type="text" value="0"/>	<input type="text" value=""/>
ISR diff%	<input type="text" value="1"/>	<input type="text" value=""/>

Quick access for number of decimal digits: right mouse button in the respective header row in VIEW.

---

**SUB Format...**

It is possible to switch between the 2 possible SUB sample display formats: SUBxxxxxxx or SUB###DG###TP###..

**Display Time as...**

Times (for SUB and VAL samples) are given either in hours or in the \_d\_h\_m or in „hh:mm“ format (50.5h = 2d2h30m = 50:30) format.

**Display LOQ as...**

<LOQ	'<LOQ' as text (or another text defined in the Reporting Texts dialog)
<LOQ (NUMERIC)	Instead of <LOQ the results show the LOQ as value taking the respective dilution factor into account. Example: if LOQ is 0.5 for DF=1: <0.5, or for DF=2 correspondingly <1.0, or for DF=0.5: <0.25 etc.
VALUE	the found concentration is displayed

**Display Not Correct SUB results as...**

There are three ways to display (yet) not correct SUB results:

:

- Result is in red color only
- Result is in red color with added **X** character
- Result is displayed as text „n.correct“ or „n.report“



### Acceptance Criteria...

The following figure shows the default values for the acceptance criteria in the individual batches. If other acceptance criteria are to be applied to the current project, they must be entered here.

OR

After a change, all batches of the current project are checked for compliance with the new acceptance criteria. If one or more analytes do not meet the new acceptance criteria the corresponding batch number is displayed with an exclamation mark '!' after the number. According to the acceptance criteria check results the respective batch tags are displayed (in the batch view) green or red. Further details on the acceptance checks are logged in the Audit Trail.

dbLabCal is able to set the E chromatogram flag and/or the S chromatogram flags automatically according to the respective users settings. This function is an extension of the acceptance criteria functionality.

User chooses the percent value of the accepted deviation. DBLABCAL checks all unknown samples (SUB) and set the flag automatically.

All E or S flags in the project can be set back to Y with the „reset!“ option. If the option is set to “off” only, current E/S flag will be NOT reset.

Acceptance Criteria		Standards	QC's
Max. deviation from nom. conc. at LOQ / LQC		20 %	15 %
Max. deviation from nom. conc. (remaining CALs/QCs)		15 %	15 %
Min. number of used and accepted values		6	4
			1
Set flag to # (all samples except CAL)if CV >		20 %	
Set flag for QCs/VALs to E if deviation higher than		reset!	?
Set flag for SUBs to S if IS not within mean(CAL) ±		reset!	S
Set flag for QCS to S if IS not within mean(CAL) ±		off	S
? At least 1 QC's at the same concentration must fit acceptance criterion!			
?? All A flags will be set automatically to J!			
???? All S flags will be set automatically to J!			
		Cancel	Ok

*Quick access: double-click on the acceptance status display in batch area.*


**Reporting Texts...**

The following dialog is available for adapting the output of the various texts to the current project.

Reporting Texts		
<b>General Texts:</b>	<b>Results:</b>	<b>Validation Samples</b>
Sequence batch	Sample destroyed during sampl.preparation destroyed	Default label Used label
Subject sub	No correct result available n.correct	NR: Matrix/Room Temperature NR: Matrix/Room Temperature
Period per	No sample NOS	NK: Matrix/Refrigerator (5±3°C) NK: Matrix/Refrigerator (5±3°C)
Time time	Not analyzed NOA	NG: Matrix/Freezer (-20±5°C) NG: Matrix/Freezer (-20±5°C)
Dilution Factor dil.f	Nothing is reported NOR	NT: Matrix/Deep Freezer (-77±5°C) NT: Matrix/Deep Freezer (-77±5°C)
	Sequence not accepted batch failed	ER: Extracts/Room Temperature ER: Extracts/Room Temperature
	Insufficient sample amount for re-assay insuff.s.	EK: Extracts/Refrigerator (5±3°C) EK: Extracts/Refrigerator (5±3°C)
	Conc. Below the Quantification Limit <LOQ	EG: Extracts/Freezer (-20±5°C) EG: Extracts/Freezer (-20±5°C)
		ET: Extracts/Deep Freezer (-77±5°C) ET: Extracts/Deep Freezer (-77±5°C)
<b>Statistics:</b>		PR: Validation Samples PR: Validation Samples
Standard deviation sd		PK: Validation Samples PK: Validation Samples
Coefficient of variation cv [%]		PG: Pools (freeze/thaw) PG: Pools (freeze/thaw)
Range range		PT: Validation Samples PT: Validation Samples
Deviation of the means bias [%]		AR: Other Matrix AR: Other Matrix
Deviation of the results dev [%]		AK: Validation Samples AK: Validation Samples
Rel.deviation from 0-Value (Stability) rel [%]		AG: Validation Samples AG: Validation Samples
<b>Explanation for CAL, QC and VAL flags:</b>		AT: Validation Samples AT: Validation Samples
N N: not used for statistical evaluation due to chromatographic error		BR: Validation Samples BR: Validation Samples
S S: not used for statistical evaluation due to not accepted IS		BK: Validation Samples BK: Validation Samples
E A: not used for statistical evaluation due to deviation from nominal conc > _1_%		BG: Validation Samples BG: Validation Samples
D D: destroyed during sample preparation		BT: Validation Samples BT: Validation Samples
X -: not measured		
		Cancel Ok

The placeholder *\_1\_* is replaced by the valid value upon printout/export.

It is the value chosen for the automatic E flag setting in the acceptance criteria dialog. If it is set to „off“, the chosen values for LOQ and lowest QC are used (see previous page).

**Set Target Folders...**

DBLABCAL expects a folder structure like:

*example*

folder	i:\unit01\	content (example)
Project01	XX001	
batch01	SP01	Chromatograms of 1st batch in plasma
batch02	SP02	Chromatograms of 2nd batch in plasma
batch03	SU01	Chromatograms of 1st batch in urine
etc...		
Project02	XX002	
batch01	AB01	
batch02	AB02	
etc...		
Project03	XX003	
etc...		

and tries to „jump“ into the folder e.g. during the import accordingly.


*This folder structure is not required by DBLABCAL. Well organized folder structure makes the daily work much easier and is strongly recommended anyway.*

**Standard Column Width!**

The column widths in the tables can be changed if required (with the mouse). By clicking on this menu option the original column widths are restored.

## Analyte

This menu option is only offered when more than one analyte was measured in a project. The project data of the selected analyte are loaded.



*Quick access: click with right mouse button on the name of the analyte.*

## Regression Model, Weighting, Readings

These menu options are accessible only if the project is a validation. After any changes to readings, weighting or regression model the data for the whole project are recalculated based on the new settings.

**These changes automatically affect all analytes that were measured in the same chromatogram! If it is necessary to change the parameter for just one analyte, it has to be done via PROJECT | EDIT...**

It is possible to change only one parameter each here. In case one has to change more/all regression parameters at once, it is more reasonable to do it in the dialog PROJECT | EDIT... See page 20. The whole project is re-calculated only once in such case

*Quick access: click with right mouse button on the corresponding terms.*

*In order to make changes to regression model, weighting or readings in a routine project the project type must first be changed from 'routine' to 'validation'. (via PROJECT | Edit... or by clicking with the right mouse button on the project type name)*

## Audit Trail and Electronic Signature

The following user actions are recorded in the audit trail with date, time and name.  
In some cases, the electronic signature is required and is also recorded in the audit trail.

Action	Electronic Signature (ES)
Defining a new project	
Editing project	X
Editing project permissions	X
Changing regression model, weighting, readings	
Changing acceptance criteria	X
Importing batch data	
Recalculating batch Checking acceptance criteria	
Changing batch status	X only if batch status changed from or to „excluded“
Changing unit or batch number	X
Changing batch start date	X
Changing sample names	X only if meanwhile batch status set
Changing a chromatogram flag	
Selecting final results of re-assays	
Changing project status	X only if project status reset to “started”
QC start and end QA end	X (when finishing QC)
Electronic signature entry	N/A

## Print...

Project results can be printed out at any time. It is recommended to use only printouts from released studies for project documentation purposes.

If several analytes (in a chromatogram) were measured in the project the user has the choice to print the selected data only for the current analyte or simultaneously for all analytes (from a chromatogram).

The following points are common to all analytes if they were measured in a chromatogram:

- BATCH LIST
- SUBJECTS (ONLY IN STANDARD FORMAT, NOT IN LIST OF DATA FORMAT)
- SAMPLES TO BE RE-ANALYZED
- RE-ASSAYED SAMPLES
- INCURRED SAMPLES
- NOT USED CALS, QCs AND VALS
- NOT USED CHROMATOGRAMS
- NOT USED BATCHS
- DATE OF ANALYSES
- AUDIT TRAIL

**Print**

Printer: \\vaidepoly\IT\_LJ4200\_S (Ne04:)

Subject Samples Results

Standard from: [dropdown]

List of Data to: [dropdown]

existing only

Statistics: Calibration Curves

Statistics: Regression Parameters

Statistics: Quality Control Samples

Samples to be Re-Analysed

Re-Assayed Samples

Incurred Samples

Not Used CALs, QCs and VALs

Not Used Sequences

Not Used Chromatograms

Date of Analyses

Hints

Audit Trail

Table of Contents/Memo

Validation Samples

NR: Matrix/Room Temperature

NK: Matrix/Refrigerator (5±3°C)

NG: Matrix/Freezer (-20±5°C)

NT: Matrix/Deep Freezer (-75±15°C)

ER: Extracts/Room Temperature

EK: Extracts/Refrigerator (5±3°C)

EG: Extracts/Freezer (-20±5°C)

ET: Extracts/Deep Freezer (-75±15°C)

PR: Validation Samples

PK: Validation Samples

PG: Pools (freeze/thaw)

PT: Validation Samples

AR: Other Matrix

AK: Validation Samples

AG: Validation Samples

AT: Validation Samples

BR: Validation Samples

BK: Validation Samples

BG: Validation Samples

BT: Validation Samples

QCs, DF<>1

Select All  print log

Print/Export Options...

Cancel Ok

Subjects are printed in the selected display format (Standard, List of Data, all or only measured samples).

Choosing "Standard" as display format allows for printing the results for the individual subjects separately by entering the number of the first and the last subject.

Choose all Options with SELECT ALL. With Keeping CTRL key pressed while clicking on SELECT ALL button will deselect all options.

The appearance of the printout, e.g. font, font size, margins, shading etc., can be modified here in the Print dialog or via the menu PROJECT | OPTIONS | PRINT/EXPORT OPTIONS.

The current standard Windows printer is displayed and can be changed here.

Audit Trail details column is occasionally shortened to fit the page width. Use export to Word or Excel (windows clipboard) to print the whole details content.

Selecting PRINT LOG will finally print a summary of the whole print job.



**Export...**

A dialog very similar to the Print dialog is displayed. MS Word does not need to be started first since DBLABCAL starts Word automatically if required.

**Export**

Subject Samples Results

Standard from: [dropdown]

List of Data to: [dropdown]

existing only

Statistics: Calibration Curves

Statistics: Regression Parameters

Statistics: Quality Control Samples

Samples to be Re-Analysed

Re-Assayed Samples

Incurred Samples (Random Repeats)

Not Used CALs and QCs

Not Used Sequences

Not Used Chromatograms

Date of Analyses

Hints

Audit Trail

Table of Contents/Memo

Validation Samples

NR: Plasma/Room Temperature

NK: Matrix/Refrigerator (5±3°C)

NG: Matrix/Freezer (-20±5°C)

NT: Matrix/Deep Freezer (-75±15°C)

ER: Prepared Samples/Room Temperature

EK: Prepared Samples/Refrigerator (5±3°C...)

EG: Prepared Samples/Autosampler(cooled)

ET: Extracts/Deep Freezer (-75±15°C)

PR: Validation Samples

PK: Validation Samples

PG: Pools (freeze/thaw)

PT: Validation Samples

AR: Matrix Effects

AK: Carry Over

AG: Validation Samples

AT: Sample Collection (Blood)

BR: Validation Samples

BK: Solutions

BG: Solutions (IS)

BT: 0h

QCs, DF<>1

Select All

Print/Export Options...

formatted

Cancel Ok

The data can be exported in a specified format, i.e. the tables are generated with margins, shading and correct column widths. The results should be exported into Word at the end of the document or, even better, into a separate document since font and font size for the Word document between the current cursor position and the end of the document may change during the export process according to the settings in DBLABCAL.

The format of the exported data, if a formatted export is desired, such as font, font size, margins, column widths, shading etc. can be modified here in the Export dialog or via the menu PROJECT | OPTIONS | PRINT/EXPORT OPTIONS.

Subjects are exported in the selected display format (Standard, Standard only measured Samples or List of Data).

*It is also possible at any time to export all currently displayed data or only the marked data via the Windows clipboard (keys CTRL-C or CTRL-INSERT, or by clicking with the right mouse button for marked data) into other programs.*

## Changing Project Status

A project, or more precisely an analyte in a project, can have one of the following statuses:

Status	Comment	Is set...
<b>Defined</b>	project was created by the analyst or PI/study director, no data imported yet	automatically on creation of project
<b>Started</b>	Data have been imported for at least one analyte	automatically on importing the 1st batch
<b>Finished</b>	All chromatographic measurements are finished	by analyst (or PI/study director)
<b>Data Released</b>	All chromatograms and batches have been accepted and results to be reported are selected - all project data locked	by PI/study director
<b>Archived</b>	Similar to 'released', however, the projects cannot be reset (un-locked) by PI/study director	by administrator

By setting the project status to '**finished**' the **analyst** confirms that all batches of the project were imported and all chromatograms of the project (with the corresponding flag) have been accepted. The PI/study director can release and lock the project data only after the analyst has declared the project as 'finished'.

By setting the project status to '**released**' the **PI/study director** locks the database. The project data can then no longer be modified. Project results are once more checked during the release process. If, for instance, the user has not yet selected final values from re-assays, the release process will be aborted.

***Also all batches are locked automatically!***

Eventually, the DB administrator sets the project status to '**archived**' (under responsibility of the archivist). Read only access to data is still allowed for all in project involved user.

**Changing the project status is an action with very far reaching consequences!**

Documentation of the QC/QA review starts after the project was released (project status = '**released**'). There are following QC/QA status' of a released project:

Status	Comment	Set by...
...	QC not started yet	
<b>QC in progress...</b>	QC started by one of the authorized QC person	Automatically after first time saving the QC comments
<b>QC-ed</b>	QC process finished (ES)	By a QC authorized user
<b>QA-ed</b>	(ES) in dbLabCal just for documentation purposes. QA details are in the QM dept.	By user with „Reviewer“ authorization (by default a user from the QM department)

### **Results for HoLaRo... / Results for Bl...**

This menu option is only displayed if unknown samples (SUB) were measured in the project, the project is released and the option has been enabled by the administrator. This option writes the results into a special ASCII file. First a dialog appears in which the data to the exported can be edited. A more detailed description is found in Appendix 1.

### **Results in ASCII File!**

This menu option is only displayed if unknown samples (SUB) were measured in the project and the project is released. This option writes the results into (an) ASCII file(s) which, for instance, can be used for the biometric evaluations. The program names and assigns the files to the analytes, displays them after the export process and enters an audit trail entry. The file name is composed of project number, name of analyte and name of matrix (examples: HA107\_ANALYT8\_PLASMA.BIO, HA107\_ANALYT8\_PLASMA-02.BIO , RT001\_SUBSTANCEX\_SERUM.BIO, etc.). If required, an automatically incremented number is added to the file name if file with the same name already exists.

The files are written into a directory set as described on page 28.

The ASCII file format is described in the dbLabCal Administration manual.

### **Exit**

Selecting Exit closes the program. The user is asked whether the current data should be loaded automatically at the next program start.

## Results Menu

Displays the results of the current project for the selected analyte and, if available, messages from the database, e.g. with regard to data discrepancies, threshold value violations in calibration, QC and/or validation samples, missing selections for values to be reported from repeat measurements etc...

The screenshot shows the dbLabCal V3 - NX023 / DAD / Seq 3 interface. The 'Results' menu is open, showing the following options:

- Statistics: Calibration Curves
- Statistics: Regression Parameters
- Statistics: Quality Control Samples
- Subject Samples Results
  - Validation Samples (selected)
  - Samples to be Re-Analysed
  - Re-Assayed Samples
  - Incurred Samples
- Not Used CALs, QCs and VALs
- Not Used Sequences
- Not Used Chromatograms
- Date of Analyses
- Hints

The 'Validation Samples' sub-menu is open, showing the following options:

- NR: Matrix/Room Temperature
- NK: Matrix/Refrigerator (5±3°C)
- NG: Matrix/Freezer (-20±5°C)
- NT: Matrix/Deep Freezer (-75±15°C)
- ER: Extracts/Room Temperature
- EK: Extracts/Refrigerator (5±3°C)
- EG: Extracts/Freezer (-20±5°C)
- ET: Extracts/Deep Freezer (-75±15°C)
- PR: Validation Samples
- PK: Validation Samples
- PG: Pools (freeze/thaw)
- PT: Validation Samples
- AR: Other Matrix
- AK: Validation Samples
- AG: Validation Samples
- AT: Validation Samples
- BR: Validation Samples
- BK: Validation Samples
- BG: Validation Samples
- BT: Validation Samples
- QCs, DF<>1

The background shows a table with columns for 'cv [%]' and 'bias [%]' and rows numbered 1 to 23.

**When loading a project or an analyte attention should first be paid to the HINTS if generated by the database.**

*For technical reasons hints to validation samples can only be displayed after the user has reviewed the data at least once.*

---

**Statistics**

... show the results of the statistical evaluations of calibration and quality control samples.

The following parameters are calculated: mean, standard deviation (sd), coefficient of variation (cv) and bias in percent.

$$mean = \frac{1}{n} \sum_{i=1}^n x_i$$


$$sd = \sqrt{\frac{\sum_{i=1}^n (x_i - \bar{x})^2}{n-1}} \quad ; n > 2 \qquad range = |x_1 - x_2| \quad ; n = 2$$

$$cv[\%] = 100 \frac{sd}{mean} \quad ; n > 2 \qquad cv[\%] = n/a \text{ (range/mean)} \quad ; n = 2$$

$$bias[\%] = 100 \frac{calc.conc - nom.conc}{nom.conc}$$

If the program finds one (and only one) batch in which the number of QC samples is higher than in all other batches, the results of the 'intra-day' statistics for that batch are calculated and displayed in addition.

If there are some QC samples with E flag (excluded from statistics by user), dbLabCal automatically adds also the statistics including those values.

The user can take individual batches for the statistical calculations or exclude them from the statistics by clicking on the button . Double-clicking or using buttons Y and N in the 'apply' columns toggles between Yes and No.

Using batches for statistics is independent of the batch status.

Only ACCEPTED batches are used for statistics by default. But it is possible to include also NOT ACCEPTED batches into the statistical evaluations.

Unit	Sequence	Comment	apply
15	1	LOQ; nr24; pg1	YES
15	2	stability stock solutions 5-ASA/Ac-5-ASA	YES
15	3	nr0; er0; ng25d	YES
15	4	stabilbity stock solution and working solution IS 4-ASA; er72; ek72	YES
15	5	intraday; recovery; pg3	YES
15	6	selectivity: six different blanks+IS/LOQ; dilution of QC 45 µg/ml (DF!	YES

Statistics:  intra day, if found     intra day, within batch  
 Accuracy instead of Bias     Accuracy instead of Dev

Show and print deviations from nominal conc.  
 Show and print calculated conc. of not used CALs/QCs too

Show and print explanation for not used results (flag: N, S, A)  
 Show and print explanation for not reported values (flag: X, D)  
 Print Not used sequences

Cancel    Ok

The user can hide parts of the statistical evaluation of the QC results and redisplay them or can initiate a printout of the respective chromatogram flag and the legend for the missing values (SHOW AND PRINT EXPLANATION FOR NOT...).


If SHOW AND PRINT DEVIATIONS... is selected the statistics obtained for calibration and QC samples include not only each individual result but also the corresponding percentage deviation from the nominal value.

Also, if SHOW AND PRINT CALCULATED.... is selected, the concentration for the Cals/QCs is displayed and printed with flags A or X in brackets. Of course, there is no calculated concentration for Cals/QCs with flags N or D.

If the project was measured at several units the user can also set the order of sorting (first by batch and then by unit number, or vice versa).

## Subject Samples Results

The results from subjects can be displayed in four different ways (Standard, Standard/measured only, List of Data, List of Data/measured only). In the standard displays the results can also be displayed graphically. The axes are scaled automatically or the user selects fixed limits for the current analyte via the menu OPTIONS | GRAPHICS... or by double-click on the graphic.



The figure above shows the data from the marked rows. A subject should be selected by clicking on the upper left corner of the subject table. A list with all subjects in the project opens. After selection of a subject (or "all" as the last entry) all results for the subject(s) are marked and the graphic is displayed. Then you can page through the subjects by using the keys PAGE↑ or PAGE↓, or through the subjects periods with CTRL - PAGE↑ or CTRL-PAGE↓.

The different dot colors have the following meaning:

- blue value Ok
- purple value Ok but should be repeated (R/C flag)
- green value is the result of repeat measurements
- red value is incorrect (further information in SAMPLES TO BE RE-ANALYZED and/or HINTS)



---

The individual periods are also identified by different line colors: black = 1st period, blue = 2nd period, green = 3rd period, magenta = 4th period, etc. The color assignments to the individual periods are shown under the graphic.

*Pressing 1-9 keys (numeric pad) changes the point size and line thickness. Keeping CTRL key simultaneously with the numeric key changes also lines colors.*

The concentration vs. time graph allows quick plausibility check of the results.


The number of expected and measured examples also helps in determining the correctness of the results obtained. (The program calculates the number of expected samples from all subjects, periods, and sampling times known so far.) If the analyst cannot explain difference between expected and measured samples it was likely due to incorrect naming of one or several subject samples. They can be found fairly quickly by "paging" through the subjects as described above.


Missing samples and samples to which the analyst cannot report a value (or does not want to) are, if present, summarized and displayed in a small table. Subject samples with an 'X' flag are also shown as 'missing' if they were only measured once.

After double-click in the 'comment' column the user can enter a comment to each measured sample. This could include, for instance, additional information to the samples, such as coding of the sponsors.

In addition to the comment written by the analyst the program always shows the so-called 'internal comment' as well. The figure indicates, for instance: first sample was OK, the sample was measured twice and the average of both measurements was used as value to be reported. The user can accept or delete the 'internal comment'. (It will nevertheless be displayed when the comment is opened again).

to be reassayed (not reliable) (mean 1+2 from 2)

In addition to each subject result the dilution factor used and the batch number(s) in which the samples were measured can be displayed in the standard display via 



### **Dilution Factor**

The dilution factor used is displayed in an additional column. Similarly, for repeatedly analyzed samples if all 're-assays' were measured with the same dilution factor. No information is available if the 're-assays' were measured with different dilution factors.

### **Batches of used chromatograms only**

The batch number(s) in which the samples used for the final results were measured are shown in an additional column.

### **Batches of all chromatograms**

Batch number(s) in which the unknown sample was measured are shown in an additional column even if the results were not used for the final result.

### **Comment**

Comment column is shown or hidden.

### **Show LOQ and 10xLOQ**

The results graphic display is shown with auxiliary lines for the quantification limit and the 10-fold value of the quantification limit.

The status line of the 'List of Data' display shows additional information to the sample under the table cursor.

## RandRep

After click on RandRep button, the RandomRepeats dialog opens. Here it is possible to choose „Random Repeats“ or “Incurred Samples” from all until now analyzed AND selected subject samples. In addition, 10 „back-up“ samples are selected. They should be used in case that it was not possible to analyze the originally selected samples e.g. due to low sample volume.

This list can be printed.

Random Repeats					
	Reserve	sub	per	time	Comment
1		1	1	2.50	
2		1	1	3.00	
3	R9	1	2	0.50	
4		1	2	2.00	
5		1	2	2.50	
6	R2	1	2	24.00	
7		2	1	0.50	
8		2	1	7.00	
9		2	2	0.50	
10		3	1	0.00	
11		3	1	3.50	
12		4	1	16.00	
13		4	2	0.00	
14		4	2	1.00	
15		4	2	2.50	
16		4	2	3.00	
17		4	2	10.00	
18		6	1	3.00	
19		6	2	3.50	
20	R1	6	2	6.00	
21	R10	7	1	2.00	
22		7	2	0.50	
23		7	2	2.00	
24		7	2	16.00	
25		7	2	240.00	
26		8	1	2.00	
27	R7	8	1	8.00	
28		8	2	2.00	
29		8	2	12.00	
30		9	1	0.50	
31		9	1	1.00	
32		9	1	3.00	
33	R6	9	1	7.00	
34		9	1	12.00	

Print    \\Aaidepawel\IT\_LJ1300N\_02\_S (Ne03:)    10 % (97 + 10)    Close

All „Random Repeats“ are marked yellow in RESULTS SUBJECT SAMPLE RESULTS, RE-ASSAYED SAMPLES and in the batch list.

## Validation Samples

There are 20 “screens” available for the validation samples. These screens may be used freely for (stability) evaluations. The screens are defined by the character combinations NR, NK,...EK...PT...AG...BG, BT. Nominal concentration of the VAL sample determines the column and the time value determines the group of samples used for statistical calculations (mean, sd, cv).


The 0 values are assigned to the individual types of stability samples as follows.

- If only one set of stability 0 values exists, it will automatically be assigned to all available “matrix/temperature” combinations.
- If one or several sets of stability 0 values with the same combination of “matrix/temperature” exist for a given combination of “matrix/temperature”, they will be displayed together.

The screenshot shows the 'Validation Samples AR: Other Matrix' window in dbLabCal V3. The window title is 'dbLabCal V3 - NX002 / Diazepam / Seq 1'. The project name is 'NX002' and the analyte is 'Diazepam / Plasma'. The validation method is 'HPLC-UV'. The table below shows the data for 38 rows, grouped by time (0h, 30m) and sequence (4, 6).

Time	Sequence	600	1100
1	0h	4	1060
2			1080
3			1090
4			1100
5	mean		1080
6	sd		13.5
7	cv [%]		1.8
8	bias [%]		-1.5
9	rel [%]		100.0
10	30m	4	1100
11			1130
12			1110
13			1130
14	mean		1120
15	sd		16.7
16	cv [%]		1.5
17	bias [%]		1.5
18	rel [%]		103.1
19	0h	6	558
20			558
21			535
22			560
23	mean		553
24	sd		12.0
25	cv [%]		2.2
26	bias [%]		-7.9
27	rel [%]		100.0
28	30m	6	538
29			537
30			533
31			534
32	mean		536
33	sd		2.65
34	cv [%]		0.5
35	bias [%]		-10.7
36	rel [%]		96.9
37			
38			

The “next” set of stability 0 values will automatically be set to 100% and the relative deviation of the stability values will be calculated in percent, based on this stability 0 set.

If no corresponding set of stability 0 values is available for a combination of “matrix/temperature” the analyst can manually assign the 0 values via the  button.

Options Validation Samples					
Default label	Used label	<input type="checkbox"/> Use stability 0 set of...		calculate with	
NR: Matrix/Room Temperature	NR: Matrix (RT) protected from light	<input type="checkbox"/>	NR: Matrix (RT) protected from light (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
NK: Matrix/Room Temperature	NK: Matrix/Room Temperature (stabilized)	<input type="checkbox"/>			<input type="checkbox"/>
NG: Matrix/Freezer (-20±5°C)	NG: Matrix/Freezer (-20±5°C)	<input type="checkbox"/>			<input type="checkbox"/>
NT: Matrix/Deep Freezer (-75±15°C)	NT: Matrix (RT) unprotected from light	<input type="checkbox"/>	NR: Matrix (RT) protected from light (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
ER: Extracts/Room Temperature	ER: Extracts (RT) protected from light	<input type="checkbox"/>	NR: Matrix (RT) protected from light (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
EK: Extracts/Refrigerator (5±3°C)	EK: Extracts/Refrigerator (5±3°C)	<input type="checkbox"/>			<input type="checkbox"/>
EG: Extracts/Autosampler (ca.10°C)	EG: Extracts/Autosampler (ca.10°C)	<input type="checkbox"/>			<input type="checkbox"/>
ET: Extracts/other conditions...light	ET: Extracts (RT) unprotected from light	<input type="checkbox"/>	NR: Matrix (RT) protected from light (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
PR: Selectivity/Matrix test (use area)	PR: Selectivity/Matrix test (use area)	<input type="checkbox"/>			<input type="checkbox"/>
PK: VAL-Samples free	PK: VAL-Samples free	<input type="checkbox"/>			<input type="checkbox"/>
PG: Freeze/Thaw	PG: Freeze/Thaw	<input type="checkbox"/>	NR: Matrix (RT) protected from light (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
PT: VAL samples free	PT: VAL samples free	<input type="checkbox"/>			<input type="checkbox"/>
AR: Whole Blood (during sample collection)	AR: Whole Blood (during sample collection)	<input type="checkbox"/>	AK: Stock Solutions Analyte (use area) (Seq. 5)	area	<input checked="" type="checkbox"/>
AK: Stock Solutions Analyte (use area)	AK: Stock Solutions Analyte (use area)	<input type="checkbox"/>	AK: Stock Solutions Analyte (use area) (Seq. 5)	area	<input checked="" type="checkbox"/>
AG: Working Solutions Analyte (use area)	AG: Working Solutions Analyte (use area)	<input type="checkbox"/>	AG: Working Solutions Analyte (use area) (Seq. 5)	area	<input checked="" type="checkbox"/>
AT: VAL samples free (0-er)	AT: VAL samples free (0-er)	<input type="checkbox"/>	AK: Stock Solutions Analyte (use area) (Seq. 5)	area	<input checked="" type="checkbox"/>
BR: Carry-over (use area)	BR: Carry-over (use area)	<input type="checkbox"/>		area	<input checked="" type="checkbox"/>
BK: Stock solution IS (use area)	BK: Stock solution IS (use area)	<input type="checkbox"/>	BK: Stock solution IS (use area) (Seq. 5)	area	<input checked="" type="checkbox"/>
BG: Working solution IS (use area)	BG: Working solution IS (use area)	<input type="checkbox"/>			<input type="checkbox"/>
BT: VAL samples free (0-er)	BT: VAL samples free (0-er)	<input type="checkbox"/>	BK: Stock solution IS (use area) (Seq. 5)	area	<input checked="" type="checkbox"/>

Show and print deviations from nominal conc.  
 use deviation from nominal concentration only

calc. conc.  
readings  
area  
height

Cancel

The description of the validation sample groups/types (used label) can be changed to better describe the stability data for other stability conditions and/or matrices than default. (e.g. „methanolic analyte solution at 5°C“, „EDTA plasma at -80°C“, „plasma with HCl at -25°C“, etc.).

Default label will be automatically used if „used label“ is deleted

The user can also specify whether deviations from the nominal value will be displayed and/or the evaluation of the stability data will be done directly via the readings used for the analyte (p.h.r. or p.a.r.), or with peak height or peak area.

If peak height or peak area is selected, dbLabCal shows automatically also internal standard results if available and possible.

Validation Samples, special case „Matrix Test“

Following example shows matrix effect evaluation for 6 different matrix samples. Each matrix was evaluated by analysing a blank sample (2 injections) and matrix sample spiked with LOQ concentration (6 injections)

Step1: define the column position in the validation result evaluation table by adding the column position using the „hidden“ position of the concentration value.

memo			
Sequence List			
	file/sampleID	name	analy
33	id_00017	w.b.	Perin
35	id_00018	val pr 1 0.10001	Perin
37	id_00019	val pr 1 0.10002	Perin
39	id_00020	val pr 1 0.10003	Perin
41		10004	Perin
43		10005	Perin
45	id_00021	val pr 1 0.10006	<b>BLANKS</b>
47	id_00024	val pr 1 0.10001	
49	id_00025	val pr 1 0.10002	
51	id_00026	val pr 1 0.10003	
53	id_00027	val pr 1 0.10004	
55	id_00028	val pr 1 0.10005	
57	id_00029	val pr 1 0.10006	Perin
59	id_00030	val pr 0 0.10001	Perin
61	id_00031	val pr 0 0.10001	Perin
63	id_00032	val pr 0 0.10001	Perin
65	id_00033	val pr 0 0.10001	Perin
67		6 matrices, 6x each 0001	Perin
69		0001	Perin
71	id_00036	val pr 0 0.10002	Perin
73	id_00037	val pr 0 0.10002	Perin
75	id_00038	val pr 0 0.10002	Perin
77	id_00039	val pr 0 0.10002	Perin
79	id_00040	val pr 0 0.10002	Perin
81	id_00041	val pr 0 0.10002	Perin
83	id_00042	val pr 0 0.10003	Perin
85	id_00043	val pr 0 0.10003	Perin
87	id_00044	val pr 0 0.10003	Perin
89	id_00045	val pr 0 0.10003	Perin
91	id_00046	val pr 0 0.10003	Perin
93	id_00047	val pr 0 0.10003	Perin
95	id_00048	val pr 0 0.10004	<b>LOQS</b>
97	id_00049	val pr 0 0.10004	
99	id_00050	val pr 0 0.10004	
101	id_00051	val pr 0 0.10004	
103	id_00052	val pr 0 0.10004	
105	id_00053	val pr 0 0.10004	
107	id_00054	val pr 0 0.10005	Perin
109	id_00055	val pr 0 0.10005	Perin
111	id_00056	val pr 0 0.10005	Perin
113	id_00057	val pr 0 0.10005	Perin
115	id_00058	val pr 0 0.10005	Perin
117	id_00059	val pr 0 0.10005	Perin
119	id_00060	val pr 0 0.10006	Perin
121	id_00061	val pr 0 0.10006	Perin
123	id_00062	val pr 0 0.10006	Perin
125	id_00063	val pr 0 0.10006	Perin
127	id_00064	val pr 0 0.10006	Perin
129	id_00065	val pr 0 0.10006	Perin
131	id_00066	w.b.	Perin

Step 2: enter the text „matrix“ in „Used label“ field and set „calculate with“ to „area“ in Option Validation Samples dialog.

**Options Validation Samples**

Default label	Used label	Use stability 0 set of...	calculate with	
NR: Matrix/Room Temperature	NR: Matrix/Room Temperature			<input type="checkbox"/>
NK: Matrix/Refrigerator (5±3°C)	NK: Matrix/Refrigerator (5±3°C)			<input type="checkbox"/>
NG: Matrix/Freezer (-20±5°C)	NG: Matrix/Freezer (-20±5°C)			<input type="checkbox"/>
NT: Matrix/Deep Freezer (-77±5°C)	NT: Matrix/Deep Freezer (-77±5°C)			<input type="checkbox"/>
ER: Extracts/Room Temperature	ER: Extracts/Autosampler (10°C)	ER: Extracts/Autosampler (10°C) (Seq. 3)	area	<input checked="" type="checkbox"/>
EK: Extracts/Refrigerator (5±3°C)	EK: Extracts/Refrigerator (5±3°C)			<input type="checkbox"/>
EG: Extracts/Freezer (-20±5°C)	EG: Extracts/Freezer (-20±5°C)			<input type="checkbox"/>
ET: Extracts/Deep Freezer (-77±5°C)	ET: Extracts/Deep Freezer (-77±5°C)			<input type="checkbox"/>
PR: Validation Samples	PR: Selectivity/Matrix Test	PR: Selectivity/Matrix Test (Seq. 3)	area	<input type="checkbox"/>
PK: Validation Samples	PK: Validation Samples			<input type="checkbox"/>
PG: Pools (freeze/thaw)	PG: Pools (freeze/thaw)			<input type="checkbox"/>
PT: Validation Samples	PT: Validation Samples			<input type="checkbox"/>
AR: Other Matrix	AR: Reliability	AR: Reliability (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
AK: Validation Samples	AK: Specificity_without Glucuronides	AK: Specificity_without Glucuronides(Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
AG: Validation Samples	AG: Specificity_with 100ng/mL Y130	AG: Specificity_with 100ng/mL Y130(Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
AT: Validation Samples	AT: Validation Samples			<input type="checkbox"/>
BR: Validation Samples	BR: Carry-over		area	<input checked="" type="checkbox"/>
BK: Validation Samples	BK: Validation Samples			<input type="checkbox"/>
BG: Validation Samples	BG: Validation Samples			<input type="checkbox"/>
BT: Validation Samples	BT: Validation Samples			<input type="checkbox"/>

Show and print deviations from nominal conc.  
 use deviation from nominal concentration only

Cancel Ok

**Result:**

dbLabCal V3 - QA639 / Seq. 1

Project Results Sequence Seq. Data Chromatography db QC+QA Extra

**Project**  
 QA639 sponsorcode xcgbgfhhfghv This project is part of a Multi-Site GLP Study  
 Pe and Pet in Human Plasma (pediatric samples) DEV  
 memo IT Vagaday Vagaday

**Analyte**  
 Validation

**Validation Samples PR: Selectivity/Matrix Test ( peak area )**


#	Sequence	0.100	0.100	0.100	0.100	0.100	0.100
1	LOG	12255	12276	11869	16510	13189	11392
2		13031	12127	11652	15333	12721	12248
3		12238	12125	11732	14959	13577	11234
4		12936	11331	12030	15392	12795	11223
5		12766	11781	11911	14284	14101	11777
6		12632	12027	12132	15384	12872	12228
7	mean	12643	11945	11888	15310	13209	11684
8	sd	336	342	179	725	539	474
9	cv [%]	2.66	2.87	1.51	4.73	4.08	4.06
10		-	-	-	-	-	-
11	rel [%]	100.00	100.00	100.00	100.00	100.00	100.00
12	BLANK	836	1120	553	3375	1650	997
13		432	558	397	2385	1522	766
14	mean	634	839	475	2880	1586	881
15	range	404	562	155	990	128	231
16	range [%]	63.66	67.06	32.72	34.37	8.07	26.26
17		-	-	-	-	-	-
18	rel [%]	5.01	7.02	4.00	18.81	12.01	7.54
19							
20							

## Validation Samples, special case „Carry-Over Test“

Following example shows how carry-over was evaluated in 6 different batches by injecting the samples in the sequence ULQ->BLANK->LOQ

*It is also possible to inject more than one BLANK sample e.g. ULQ->BLANK-> BLANK->LOQ...*

Step1:



dbLabCal - FX064 / ... / Seq 1

Project Results Sequence Seq. Data Chromato

**Project**  
QA639 sponsorcode xcgbgfh  
P and Pat in Human Plasma (pediatric samples) DEV  
memo IT Vagada

**Sequence**  
1 UNIT01 Vagada  
Carry-over(1)/Between-assay(1)/Dilution Test  
memo

**Sequence List**

	file/sampleID	name
1	id_00001	5ul k-ref
2		
3	id_00002	w.b.
4		
5	id_00003	val br 0h 50.0
6		val br 0h 50.0
7	<del>id_00004</del>	val br 0h 0.00
8		val br 0h 0.00
9	id_00005	val br 0h 0.100
10		val br 0h 0.300
11	id_00006	w.b.
12		



Step 2: enter the text „carry“ and „over“ in „Used label“ field and set „calculate with“ to „area“ in Option Validation Samples dialog.

**Options Validation Samples**

Default label	Used label	Use stability 0 set of...	calculate with	
NR: Matrix/Room Temperature	NR: Matrix/Room Temperature			<input type="checkbox"/>
NK: Matrix/Refrigerator (5±3°C)	NK: Matrix/Refrigerator (5±3°C)			<input type="checkbox"/>
NG: Matrix/Freezer (-20±5°C)	NG: Matrix/Freezer (-20±5°C)			<input type="checkbox"/>
NT: Matrix/Deep Freezer (-77±5°C)	NT: Matrix/Deep Freezer (-77±5°C)			<input type="checkbox"/>
ER: Extracts/Room Temperature	ER: Extracts/Autosampler (10°C)	ER: Extracts/Autosampler (10°C) (Seq. 3)	area	<input checked="" type="checkbox"/>
EK: Extracts/Refrigerator (5±3°C)	EK: Extracts/Refrigerator (5±3°C)			<input type="checkbox"/>
EG: Extracts/Freezer (-20±5°C)	EG: Extracts/Freezer (-20±5°C)			<input type="checkbox"/>
ET: Extracts/Deep Freezer (-77±5°C)	ET: Extracts/Deep Freezer (-77±5°C)			<input type="checkbox"/>
PR: Validation Samples	PR: Selectivity/Matrix Test	PR: Selectivity/Matrix Test (Seq. 3)	area	<input type="checkbox"/>
PK: Validation Samples	PK: Validation Samples			<input type="checkbox"/>
PG: Pools (freeze/thaw)	PG: Pools (freeze/thaw)			<input type="checkbox"/>
PT: Validation Samples	PT: Validation Samples			<input type="checkbox"/>
AR: Other Matrix	AR: Reliability	AR: Reliability (Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
AK: Validation Samples	AK: Specificity_without Glucuronides	AK: Specificity_without Glucuronides(Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
AG: Validation Samples	AG: Specificity_with 100ng/mL Y130	AG: Specificity_with 100ng/mL Y13(Seq. 2)	calc. conc.	<input checked="" type="checkbox"/>
AT: Validation Samples	AT: Validation Samples			<input type="checkbox"/>
BR: Validation Samples	BR: Carry-over		area	<input checked="" type="checkbox"/>
BK: Validation Samples	BK: Validation Samples			<input type="checkbox"/>
BG: Validation Samples	BG: Validation Samples			<input type="checkbox"/>
BT: Validation Samples	BT: Validation Samples			<input type="checkbox"/>

Show and print deviations from nominal conc.  
 use deviation from nominal concentration only

Cancel Ok

Result:

dbLabCal V3 - QA639 / / / Seq. 2

Project Results Sequence Seq. Data Chromatography db QC+QA Extra

**Project**  
 QA639 sponsorcode xcgbgfhbghv This project is part of a Multi-Site GLP Study  
 Pe and Pet in Human Plasma (pediatric samples) DEV  
 memo IT Vagaday Vagaday

**Validation Samples BR: Carry-over ( peak area )**

	Sequence	0.00	0.100	50.0
1	1 *	1230	11661	5615080
2	2 *	3076	6363	3452020
3	3 *	2398	10808	5449620
4	4 *	2919	12486	5173020
5	5 *	2102	13416	5841030
6	mean	2345	10947	5106154
7	sd	736	2739	956249
8	cv [%]	31.40	25.02	18.73
9	blank in % LOQ	21.42	-	-
10	blank in % UOQ	0.05	-	-
11				
12				
13				

## Re-Assays

### Samples to be Re-analyzed

SAMPLES TO BE RE-ANALYZED shows all subject samples (chromatograms) with a flag 'N', 'R'; 'C', 'D', 'V', '<', '>' or '+' which have not been re-analyzed yet. This list includes all chromatograms where the flag of any analyte in the chromatogram meets the above mentioned condition. Samples with the 'X' flag are ignored.

sub	per	time	Analyte	ng/mL (calc)	dil	Seq	Details	OK
1	1	12.00	Analyte1	38.5	1.0	4	DeptMan: chromatogr. error	
2	2	9.00	Analyte2	72.9	1.0	4	DeptMan: internal standard area/height not accepted	
3								
4								
5								
6								
7								

The column 'Details' shows the 'translation' of the respective chromatogram flag and the name of the analyst (for flags: 'N', 'R', 'C') or DBLABCAL (for flags: '<', '>', '+'). If a chromatogram flag was changed in VIEW | LIST the analyst should give a detailed reason in the 'comment' column. This comment to the chromatogram is also shown in here the 'Details' column.

The analyst can accept the sample to be analyzed shown in the last column by double-click to confirm that the sample does not need to be re-analyzed. This, of course, is only possible if DBLABCAL has set the chromatogram flag and not the analyst.

## Re-assayed Samples

RE-ASSAYED SAMPLES shows by default all standard re-assays without incurred samples.

Samples re-assayed for “standard reason” (flag N, R, D, >, <...) as well as for ISR are shown in both lists (RE-ASSAYED SAMPLES and INCURRED SAMPLES) with the respective back color (yellow or white).

Percent deviation of the re-assays is calculated and shown in the last column (sequence: 1.-2. | 2.-3. | 3.-1. measurement) and expressed in percent based on the respective mean. (See also OPTIONS RE-ASSAYS). If the individual values deviate from each other by more than specified in the dialog OPTIONS RE-ASSAYS the mean cannot be used for final result. (but...remember the study director special wish key...)

sub	per	time	1st meas.	Seq	flag	2nd meas.	Seq	flag	3rd meas.	Seq	flag	reported	Choice	Analyte	diff(%)
1	5	1	12.00	305	6	R	308	11	Y	-	-	305	M12	Analyte1	11H
2				50.7	6	R	51.0	11	Y	-	-	50.8	M12	Analyte2	0H
3	5	1	14.00	341	6	R	333	11	Y	-	-	337	M12	Analyte2	2H
4				83.1	6	R	80.6	11	Y	-	-	81.8	M12	Analyte2	3H
5	5	1	16.00	442	6	R	423	11	Y	-	-	432	M12	Analyte1	4H
6				111	6	R	104	11	Y	-	-	107		Analyte2	6H
7	5	1	24.00	242	6	R	247	11	Y	-	-	244		Analyte1	2H
8				43.2	6	R	44.2	11	Y	-	-	43.7		Analyte2	2H
9	7	1	16.00	300	7	R	283	11	Y	-	-	295		Analyte1	1H
10				43.7	7	R	40.1	11	Y	-	-	41.9		Analyte2	1H
11	7	1	24.00	522	7	R	498	11	Y	-	-	510		Analyte1	7H
12				75.6	7	R	71.3	11	Y	-	-	73.5		Analyte2	6H
13	9	1	16.00	194	8	R	172	11	Y	-	-	183		Analyte1	12H
14				48.7	8	R	42.0	11	Y	-	-	45.4		Analyte2	15H
15	9	1	24.00	371	8	R	347	11	Y	-	-	359	n.report.	Analyte1	7H
16				120	8	R	111	11	Y	-	-	116		Analyte2	8H
17	14	2	11.00	145	10	R	146	11	Y	-	-	146		Analyte1	11H
18				21.4	10	R	21.6	11	Y	-	-	21.5		Analyte2	11H
19	14	2	12.00	164	10	R	167	11	Y	-	-	166		Analyte1	2H
20				26.6	10	R	27.1	11	Y	-	-	26.9		Analyte2	2H
21	14	2	14.00	253	10	R	255	11	Y	-	-	254		Analyte1	1H

The analyst selects (right mouse button in the ‘Choice’ column) the value to be reported from re-assays. DBLABCAL shows only the allowed/accepted choices. E.g., if there is only one choice for value to be reported (e.g. 1st measurement=flag ‘N’, ‘>’, ‘<’, ‘D’ etc., 2. measurement=flag ‘J/Y’), only this one is allowed. It is therefore recommended, and it saves several clicks, to let the reportable values be entered **automatically** first (see below).

There is also the option of not reporting any value. In this case a reason should be given in the ‘comment’ column in RESULTS | SUBJECT SAMPLES.


Incurred sample results are shown in a separate list/screen. Working with Incurred samples and the standard re-assay results is identical.

After the final results from re-assay was selected, the 'Choice' column shows a code. The meanings of these codes are explained in the following table:

Short	Meaning
1	<b>The result of the first measurement is reported</b> - no correct second (or third) measurement
2a 2b 2c 2d	<b>The result of the second measurement is reported</b> - no correct first (and third) measurement available - no reliable first measurement available - first measurement above the calibrated range - first measurement below the reduced calibrated range
3	<b>The result of the third – ninth measurement is reported</b> - no correct first and second measurements available
M12 M13 M23 M123 M14 M24 M34 M124 M134 M234 M1234	<b>The mean of two/three/four measurements is reported</b> - first measurement confirmed by second measurement *) - first measurement confirmed by third measurement *) - second measurement confirmed by third measurement *) - first and second measurements confirmed by third measurement *) Etc.

\*) Difference of individual measurements is less than the value selected in the dialog OPTIONS | RE-ASSAYS.


After selecting the first submenu, or if this selection was done according to Lang&Bolton, the 'Choice' column does not show any of the above codes but only the information which result is reported (1st value, 2nd value, mean 1+2, mean123 etc.)

After click on  the Re-assay Option dialog is shown. In this dialog the PI/study director can chose the re-assay reporting rule, set the acceptance criteria for incurred samples as well as start automatic choice of final result.

Following reporting rules are possible:

- dbLabCal rule (confirmed mean, otherwise next assay is required ...)
- first accepted result is reported
- according Lang/Bolton rule  
(J.R.Lang and S.Bolton, J.Pharm.Biomed.Anal. **9**, 357-361 (1991), see Att. 14)
- according HoLaRo rule, see Att. 15

PI/study director also define the accepted deviation [%] of the individual results from their mean. The user can define various acceptable deviations. If a result is less than or equal to the lowest QC the corresponding criterion is applied. Both criteria are set to 28% by default when defining a new project.



The automatic choice of final results is controlled in the **Action** part of the dialog. Either the not yet set re-assays (2nd checkbox) can be set automatically or all „re-assays“ are set automatically (3rd checkbox). In that case dbLabCal will „overwrite“ already user set choice or even remove the currently set choice in case no distinct result exists.

If only one „valid“ result is available from the samples that were re-assayed, it will be automatically reported.

Criterion for selecting final result from two (or three) valid measurements (chromatogram flag: J, R/C) is the difference of the (two) values expressed in percent of their mean value. If the difference is less than the value selected in OPTIONS REASSAY the mean is calculated, otherwise no automatic selection occurs.

At first accepted result rule the first valid assay (Chromatogramm flag: J, O, W/R) is reported.

Flow chart of the Lang&Bolton and HoLaRo rules are in attachments 14 and 15, respectively.

If more than maximally allowed measurements (3, 4 or max.9) were found (see HINTS) the analyst or the PI/study director must decide which three or four chromatograms are to be evaluated. The other chromatograms should be renamed (e.g. to DIV).

Multiple measurements where all concentrations were <LOQ are also shown here as an overview. If all LOQs (after taking the dilution factor into account) were identical, <LOQ is automatically taken as value to be reported and cannot be changed. If the LOQs, after taking the dilution factor into account, are different the analyst must decide which LOQ value will be reported.

It is not possible to choose final result, if the cell's back color is grey, It is the case if all concentrations were <LOQ or if another user with higher authorization have made the choice. The user name who have choosen the final result is shown automatically in the tooltips

## Incurring Samples

In this list re-assays performed for ISR are shown.

Samples re-assayed for “standard reason” (flag N, R, D, >, <...) as well as for ISR are shown in both lists (RE-ASSAYED SAMPLES and INCURRED SAMPLES) with the respective back color (yellow or white).

**Options Reassays**

Rules for reported Result (Reassays)

- dbLabCal default (mean if diff <= 28% \*), total 2 or 3 assays
- 1st accepted result (total 2 to 9 assays)
- Lang & Bolton Rule (median/mw, <15% / <30%, total 2 or 3 assays)
- HoLaRo Rule (mean if diff. <= 28% \*), total 2 or 3 assays

Max. deviation between re-assays in [%] of mean: conc. > lowest QC 28 10-40%  
 conc. <= lowest QC 28 10-60%

Rules dev% calculation (Incurred Samples)

- FDA: diff. between ISR and FIRST in % of mean dev% = (ISR-A) / mean(ISR+A)
- EBF: diff. between ISR and FIRST in % of FIRST dev% = (ISR-A) / A

Max. difference between Incurred Samples in [%] 20 10-40%  
 accepted Incurred Samples at least 66.7 66.7-100%

Action: all peaks

- Reset all choices!
- Do auto choice, do not overwrite already set values!
- Do auto choice for all reassays! Overwrite already set values!
- Sort reassays according to batch number (default sorting is according batch start date)

Ok Cancel

Rules and acceptance criteria for Incurred Samples acceptance are defined in the middle part of the dialog.

Accepted deviation valid for incurred samples is in general different than the “standard re-assay” deviation (20%, 28% respectively).

Furthermore, dbLabCal checks if at least 2/3 of the Incurred Samples fulfilled the max. accepted deviation criterion

## Excluded Values

This includes

- all NOT USED CALs, QCs and VALs,  
i.e. calibration, QC or validation samples with flags 'N' or 'A/E',
- NOT USED BATCHES (status BATCH NOT ACCEPTED or EXCLUDED) and
- NOT USED CHROMATOGRAMS ('X' flag).

## Date of Analyses

The menu option DATE OF ANALYSES shows the project performance. Extraction dates, start and end of batch measurements and the comments to each batch can be edited here. In addition, the unit and batch numbers can be directly modified.


*The comment should also include information to the samples measured in the batch, e.g. 'subjects 5 and 6', Cal/QC set numbers, 'stability 48h' etc., but not the batch number.*

If the project was measured at several units the user can specify the sorting batch (first by batch number and then by unit number, or vice versa) by clicking on >>>.

DBLABCAL accept also date-related terms, such as "today, yesterday, day before yesterday, x days ago, Monday (or Mo), Tuesday (Tu)" etc. – like everywhere in the program where date-related information is required.

## Batch Menu


Batches and batch data are managed via the Batch menu.






## Import File...

A dialog appears for selecting the batch ASCII file to be imported.



User can save currently selected data folder as start folder for next imports by right mouse click on the folder list:



Currently selected data folder is saved as start folder for next imports after successful import automatically.

All available information to the selected file is shown on the right side of the dialog.

Prior to importing an ASCII file a dialog with the content of the file to be imported is shown:

Import Check											
	file/sampleID	name	analyte	rt	area	ht	rt(is)	area(is)	ht(is)	dil.f	flag
1	+	sh1001	k.ref. 20µl	Analyte1	2.41	493726	38059	2.68	535978	38989	1. Y
1		sh1001	k.ref. 20µl	Analyte2	2.68	428920	31263	2.68	535978	38989	1. Y
2	+	sh1002	w.b.	Analyte1	0	0	0	0	0	0	1. Y
2		sh1002	w.b.	Analyte2	0	0	0	0	0	0	1. Y
3	+	sh1003	blank o. i.s.	Analyte1	0	0	0	0	0	0	1. Y
3		sh1003	blank o. i.s.	Analyte2	0	0	0	0	0	0	1. Y
4	+	sh1004	blank m. i.s.	Analyte1	0	0	0	2.71	42152	2994	1. Y
4		sh1004	blank m. i.s.	Analyte2	0	0	0	2.71	42152	2994	1. Y
5	+	sh1005	cal 0.250	Analyte1	2.53	1231	101	2.8	42613	3264	1. Y
5		sh1005	cal 0.250	Analyte2	2.77	1095	83	2.8	42613	3264	1. Y
6	+	sh1006	cal 0.250	Analyte1	2.47	1182	90	2.74	42494	3109	1. Y
6		sh1006	cal 0.250	Analyte2	2.74	954	71	2.74	42494	3109	1. Y
7	+	sh1007	cal 1.00	Analyte1	2.47	4417	337	2.74	42405	3175	1. Y
7		sh1007	cal 1.00	Analyte2	2.74	3606	270	2.74	42405	3175	1. Y
8	+	sh1008	cal 1.00	Analyte1	2.5	5125	411	2.77	51159	3872	1. Y
8		sh1008	cal 1.00	Analyte2	2.77	4610	327	2.77	51159	3872	1. Y
9	+	sh1009	cal 2.50	Analyte1	2.47	13680	1044	2.74	49265	3678	1. Y
9		sh1009	cal 2.50	Analyte2	2.74	11563	858	2.74	49265	3678	1. Y
10	+	sh1010	cal 2.50	Analyte1	2.47	12969	998	2.74	48420	3641	1. Y
10		sh1010	cal 2.50	Analyte2	2.74	11140	825	2.74	48420	3641	1. Y
11	+	sh1011	cal 10.0	Analyte1	2.47	58137	4463	2.74	53476	3997	1. Y
11		sh1011	cal 10.0	Analyte2	2.77	49491	3662	2.74	53476	3997	1. Y
12	+	sh1012	cal 10.0	Analyte1	2.47	49646	3746	2.74	45241	3353	1. Y
12		sh1012	cal 10.0	Analyte2	2.74	41458	3034	2.74	45241	3353	1. Y
13	+	sh1013	cal 25.0	Analyte1	2.47	135269	10258	2.77	44775	3308	1. Y
13		sh1013	cal 25.0	Analyte2	2.77	111572	8482	2.77	44775	3308	1. Y

In this dialog, it is possible to include / exclude single chromatograms from import by clicking on + sign in the first column (include: +, exclude: -).

If more than one analyte were measured in a chromatogram double-clicking in the 'analyte' column reduces the display to only one analyte. A second double-click in the 'analyte' column redispays all analytes.

Import Check											
	file/sampleID	name	analyte	rt	area	ht	rt(is)	area(is)	ht(is)	dil.f	flag
1		sh1001	k.ref. 20µl	Analyte2	2.68	428920	31263	2.68	535978	38989	1. Y
2		sh1002	w.b.	Analyte2	0	0	0	0	0	0	1. Y
3		sh1003	blank o. i.s.	Analyte2	0	0	0	0	0	0	1. Y
4		sh1004	blank m. i.s.	Analyte2	0	0	0	2.71	42152	2994	1. Y
5		sh1005	cal 0.250	Analyte2	2.77	1095	83	2.8	42613	3264	1. Y
6		sh1006	cal 0.250	Analyte2	2.74	954	71	2.74	42494	3109	1. Y
7		sh1007	cal 1.00	Analyte2	2.74	3606	270	2.74	42405	3175	1. Y
8		sh1008	cal 1.00	Analyte2	2.77	4610	327	2.77	51159	3872	1. Y
9		sh1009	cal 2.50	Analyte2	2.74	11563	858	2.74	49265	3678	1. Y
10		sh1010	cal 2.50	Analyte2	2.74	11140	825	2.74	48420	3641	1. Y
11		sh1011	cal 10.0	Analyte2	2.77	49491	3662	2.74	53476	3997	1. Y
12		sh1012	cal 10.0	Analyte2	2.74	41458	3034	2.74	45241	3353	1. Y
13		sh1013	cal 25.0	Analyte2	2.77	111572	8482	2.77	44775	3308	1. Y
14		sh1014	cal 25.0	Analyte2	2.77	96750	7264	2.74	39214	2841	1. Y
15		sh1015	cal 100	Analyte2	2.77	422737	31247	2.77	44943	3272	1. Y
16		sh1016	cal 100	Analyte2	2.77	358738	26598	2.77	37186	2717	1. Y
17		sh1017	cal 250	Analyte2	2.77	874303	62431	2.77	39686	2838	1. Y
18		sh1018	cal 250	Analyte2	2.89	748625	50667	2.89	31794	2402	1. Y
19		sh1019	w.b.	Analyte2	0	0	0	0	0	0	1. Y
20		sh1020	qcs 0.250	Analyte2	2.8	780	58	2.8	32111	2338	1. Y
21		sh1021	qcs 0.250	Analyte2	2.77	800	58	2.77	33669	2430	1. Y
22		sh1022	qcs 0.750	Analyte2	2.8	2515	171	2.8	35030	2472	1. Y
23		sh1023	qcs 0.750	Analyte2	2.89	2280	123	2.89	32080	2175	1. Y
24		sh1024	qcs 20.0	Analyte2	2.8	58225	3907	2.8	30924	2095	1. Y
25		sh1025	qcs 20.0	Analyte2	2.92	61753	3670	2.89	32453	2189	1. Y

After click on AUTO button, DBLABCAL suggests the required batch data like batch number, unit, analyst name as well as the peaks assignments as good as possible.

The user has to check the proposed values and enter the missing data. DBLABCAL supports the user in that OK button is allowed only if the batch data are complete and reasonable.

**It is important that all files for analytes analyzed in one chromatogram be imported TOGETHER at the first batch import. This way, the information, how many and which analytes belongs together (were analyzed in one run) is saved in dbLabCal!**

DBLABCAL allows import (OK button) at next batch imports only if the analyte count in the file fits to the database known analyte count.

**It is almost impossible to import an file into an 'wrong' project.**

### Import from Empower2...

It may be possible to import data from Empower2 into dbLabCal directly, if Empower2 and dbLabCal are installed on a PC, and if the administrator has allowed this function.

The user has to log-in with his Empower2 account into Empower2 in this dialog. He loads a Project, one or more Sample Sets and finally the Results. After completion the missing data entry (button AUTO) the selected Results (+) are imported into dbLabCal.

Database: **NEUEMP2** UserID/PWD: **vagadaym / \*\*\*\*\*** login Projects: **TA348\_Acyclovir** Launch Empower

Show Sample Sets for selected Project: 1 sample set(s) selected

SampleSetName	SampleSetStartDate	SampleSetFinishDate	SystemName
run17_unit03	24-Oct-2009 06:11:33	25-Oct-2009 05:42:15	MP03
run18_unit07	23-Oct-2009 21:25:20	24-Oct-2009 23:07:44	MP07
run16_unit03	23-Oct-2009 05:40:51	24-Oct-2009 06:11:32	MP03
run15_unit07	22-Oct-2009 19:41:58	23-Oct-2009 21:25:19	MP07
run14_unit03	22-Oct-2009 07:44:45	23-Oct-2009 01:27:26	MP03
run13_unit07	21-Oct-2009 19:43:38	22-Oct-2009 19:41:57	MP07
run12_unit03	21-Oct-2009 06:40:23	22-Oct-2009 07:44:44	MP03
run11_unit07	20-Oct-2009 18:00:17	21-Oct-2009 19:43:37	MP07
run10_unit03	20-Oct-2009 06:24:10	21-Oct-2009 06:40:22	MP03

Show latest Results for selected Sample Sets: 90 chrm(s) selected

Vial	SType	Info_Sub	Period	STime	Temp	Matrx	Dilution	DateAcquired	ResultId
1	+	1 k-ref					1,0000	24-Oct-2009 06:28:27	7754
2	+	2 w.b.					1,0000	24-Oct-2009 06:44:49	7674
3	+	3 blank					1,0000	24-Oct-2009 07:01:11	7675
4	+	4 blank + IS					1,0000	24-Oct-2009 07:17:31	7676
5	+	5 cal 10.0 ng/mL					1,0000	24-Oct-2009 07:33:50	7672
6	+	6 cal 20.0 ng/mL					1,0000	24-Oct-2009 07:50:10	7664
7	+	7 cal 40.0 ng/mL					1,0000	24-Oct-2009 08:06:32	7665
8	+	8 cal 80.0 ng/mL					1,0000	24-Oct-2009 08:22:52	7666
9	+	9 cal 100 ng/mL					1,0000	24-Oct-2009 08:39:12	7667
10	+	10 cal 200 ng/mL					1,0000	24-Oct-2009 08:55:33	7668
11	+	11 cal 400 ng/mL					1,0000	24-Oct-2009 09:11:53	7669
12	+	12 cal 800 ng/mL					1,0000	24-Oct-2009 09:28:13	7670
13	+	13 cal 1000 ng/mL					1,0000	24-Oct-2009 09:44:34	7671
14	+	14 w.b.					1,0000	24-Oct-2009 10:00:54	7677
15	+	15 sub 23	1	0h			1,0000	24-Oct-2009 10:17:14	7755
16	+	16 sub 23	1	0.5h			1,0000	24-Oct-2009 10:33:33	7756
17	+	17 sub 23	1	1.0h			1,0000	24-Oct-2009 10:49:54	7757
18	+	18 sub 23	1	1.5h			1,0000	24-Oct-2009 11:06:14	7681
19	+	19 sub 23	1	2.0h			1,0000	24-Oct-2009 11:22:34	7758
20	+	20 sub 23	1	2.5h			1,0000	24-Oct-2009 11:38:54	7759
21	+	21 sub 23	1	3.0h			1,0000	24-Oct-2009 11:55:15	7760
22	+	22 sub 23	1	3.5h			1,0000	24-Oct-2009 12:11:34	7761
23	+	23 sub 23	1	4.0h			1,0000	24-Oct-2009 12:27:54	7762
24	+	24 sub 23	1	5.0h			1,0000	24-Oct-2009 12:44:14	7763
25	+	25 sub 23	1	6.0h			1,0000	24-Oct-2009 13:00:34	7764
26	+	26 sub 23	1	8.0h			1,0000	24-Oct-2009 13:16:54	7765
27	+	27 sub 23	1	10h			1,0000	24-Oct-2009 13:33:14	7766
28	+	28 sub 23	1	12h			1,0000	24-Oct-2009 13:49:37	7767
29	+	29 sub 23	1	16h			1,0000	24-Oct-2009 14:05:59	7768
30	+	30 sub 23	1	24h			1,0000	24-Oct-2009 14:22:20	7693
31	+	31 w.b.					1,0000	24-Oct-2009 14:38:40	7694
32	+	32 sub 23	2	0h			1,0000	24-Oct-2009 14:55:00	7769
33	+	33 sub 23	2	0.5h			1,0000	24-Oct-2009 15:11:20	7770


Unit: auto Project: TA348\_Acyclovir  
 MP03 HPLC MP03 Unit  
 Analyst Wick Batch No  
 Comment  
 Batch No: 17  overwrite  
 Confirm to overwrite batch! status - batch status not set yet  
 Analyst extracted started: 24-10-2009 06:28:00  
 finished: 25-10-2009 05:26:00  
 Comment to Sequence: sample set name=run17\_unit03 (id 7078)  
 ChrmCount: 89  
 Peak Assignment: Acyclovir / Plasma  
 Peak 1 Acyclovir

Cancel OK

To open the current Empower2 Project click LAUNCH EMPOWER button.

**Load...**

A list of all batches in the current project (and for the current analyte) is displayed.



Batches are sorted by unit number and batch number.

The batch number is (sometimes) followed by a character which indicates the batch status:

- \* status undetermined
  - batch accepted
  - N batch not accepted
  - X batch excluded
  - D batch marked as deleted (overwritten)
  - L batch was locked by PI
- or project was released and then project status reset to "started" since dbLabCal locks all batches when status is set to released automatically*

The meaning of batch status is described further below.

An exclamation mark (!) after the batch number indicates that (at least) one analyte did not meet the acceptance criteria applicable to the current project.


The display of overwritten and excluded batches can be activated via the option SHOW NOT USED BATCHES.

Quick access: double-clicking anywhere within the screen section 'batch'.

*The next or the previous batch can also be loaded directly via "paging" (keys CTRL + or CTRL -; PLUS/MINUS keys on the numerical key pad - in a batch view) without first opening the dialog LOAD BATCH...*

**Print...**

Data from the current batch or from all other batches can be printed irrespective of their status. However, for documentation purposes, only printouts of batches with a project status **,released'** should be used. If several analytes in a chromatogram were measured in a project the user can print the selected data either for the current analyte only or simultaneously for all analytes that were measured together with the current analyte in a chromatogram.



*Printouts from calibration and quality control samples and the complete batch list are sufficient for documentation of the batch data. Data to subject or validation samples are documented on the printouts of the project results..*

Select all options with SELECT ALL, or de-select all options with CTRL- SELECT ALL.

The appearance of the printout, e.g. font, font size, margins, shading etc., can be modified here in the Print dialog or via the menu PROJECT | OPTIONS | PRINT/EXPORT OPTIONS.

The current standard Windows printer is displayed and can be changed here as well.

**Export...**

Same procedure as described for BATCH | PRINT. Only data for the current analyte can be exported.

The screenshot shows the 'Export' dialog box with the following settings:


- Sequence:**
  - Current Sequence
  - All Sequences
  - from: 1 to: 11
  - Unit: all
  - OK Sequences Only
- Options:**
  - Calibration Curve
  - Quality Control Samples
  - Subjects
  - Validation Samples
  - Batch List
- Buttons:** Select All, Print/Export Options...,  formatted, Cancel, Ok

*Furthermore, the current data can be exported at any time via the Windows clipboard (keys CTRL-C or CTRL-INSERT or with the right mouse click (for selected data) into other programs.*

## Change Batch Status...

Batch status can be changed for each analyte individually. It means, that it is possible to accept batch data for one analyte, while rejecting the data for another analyte (due to not fulfilled acceptance criteria) of the same batch.

Batch status from/to BATCH EXCLUDED is always performed for all analytes of one chromatogram.



	Text	Comment	Is set...
*	Status undetermined!	After import each batch is automatically set to status '*'	automatically
(Y)	Batch accepted	Can be set by analyst if the batch met the <b>acceptance criteria</b> . Can (theoretically) be set by PI/study director even where the acceptance criteria have not been met.	by analyst or PI/study director
N	Batch not accepted	Is set whenever the batch does not meet the valid <b>acceptance criteria</b> . The subject samples measured in this batch are shown in SAMPLES TO BE RE-ANALYZED.	by analyst or PI/study director
X	Batch excluded	Example: "system set-up" batch or a batch that was not evaluated for other reasons... This batch is ignored in the evaluation.	by analyst or PI/study director
D	Batch deleted	A new batch (with the same batch number) was imported. This batch is ignored in the evaluation.	Automatically
L	Batch locked	Project was released and then project status reset to "started" since dbLabCal keeps batch status "locked" even if project status was reset. PI/study director have to remove batch lock explicitly if needed	Automatically or PI/study director

### Changing the batch status is an action with far reaching consequences!

*(At least for the analyst...)*

If the batch is set to status BATCH ACCEPTED, DBLABCAL checks whether the calibration curve, QC samples and  $r^2$  met the current acceptance criteria.

If the acceptance criteria were not met the batch release is aborted. Only the PI/study director has authority to set such a batch to status BATCH ACCEPTED.


*Quick access: right mouse button over status....*



## Batch Data Menu

shows several views of the current batch data.


## Standars, QCs, Subjects and Validation Samples...




The user can view the calibration curve or quality control samples of the current batch with graphics. Red dots indicate that the deviation from the nominal concentration is higher than specified in the menu PROJECT | OPTIONS | ACCEPTANCE CRITERIA.

By holding down the CTRL key and clicking on the figure the calibration curve is sequentially shown in three different ways:


x: concentration  
y: readings



x: concentration, logarithmic  
y: readings/concentration



x: concentration, logarithmic  
y: readings, logarithmic



The subject samples and validation samples measured in the current batch can also be viewed.

*There is a difference between subject samples and validation samples shown here and those from the RESULTS menu. The menu RESULTS (of the project) summarizes the values to be reported (re-assays evaluated, flag 'Y' etc...) whereas the menu VIEW (of batch data) shows only the values measured in the current batch, irrespective of re-assays or flags and samples that were measured in other batches.*

In subject view, as in RESULTS | SUBJECT SAMPLES the values can also be displayed graphically. Red colored dots indicate only that the flag is not 'Y'.

*The subject results display of a batch can be used to view the results after importing the batch, for instance, in order to immediately recognize 'not plausible results' and set the 'R' flags.*

The user can edit unit number, batch number, comment, extraction date as well as date and time of start and end of the current batch In menu RESULTS | DATE OF ANALYSES. Information to the start of the batch is very important since a subsequent validation defines the set of 0 values for the stability studies and the order of re-assays for subject results.

Furthermore, by pressing the right mouse button over the first (gray) table row in BATCH DATA view the number of decimal digits for the display of values in the respective column can be quickly changed. The number of significant digits can only be changed via the menu PROJECT | OPTIONS | NUMBER FORMAT.

## Batch List

A complete list of the chromatography data of the current batch for all analytes measured in a chromatogram is displayed in the order in which they were measured. Sample name and chromatogram flags can be modified by the analyst as long as the batch status is 'not set yet'.

dbLabCal - HA107 / AC-5-ASA / Seq 4

Project Results Sequence Seq. Data Chromatography gb QC+QA

**Project**  
 HA107 98-0435-001-L1  
 Ac-5-ASA / 5-ASA  
 memo BA Vagaday Otto

**Analyte**  
 AC-5-ASA / plasma  
 Routine HPLC-FL ng/mL  
 data released

**Model**  
 y=a+bx  
 1/x  
 p.h.t.

**Sequence**  
 4 9 Otto  
 CAL 04 - Subject: 1+ 2 / Period :1+2  
 memo Acceptance cals qcs r<sup>2</sup> a  
 batch accepted

extracted 15-Mrz-1998 a = 0.00001  
 started 16-Mrz-1998 07:47 b = 0.00984  
 finished 17-Mrz-1998 08:25  
 calculated 17-Mrz-1998 10:04 r = 0.99999

**Sequence List**

file/sample	name	analyte	rt	area	ht	rt(s)	area(s)	ht(s)	p.h.r.	dil f	ng/mL	dev [%]	flag	comment
1 ro04001a	sss	AC-5-ASA	2.77	3458351	202247	6.98	1951732	82112	2.4631	1.0	250		J	
2		5-ASA	4.98	2045394	101459				1.2356		100		J	
3 ro04002a	w.b.	AC-5-ASA	-	-	-	-	-	-	-	1.0	-		J	
4		5-ASA	-	-	-	-	-	-	-	-	-		J	
5 ro04003a	Blank	AC-5-ASA	-	-	-	6.97	1991341	84158	0.0000	1.0	-0.000523		J	
6		5-ASA	-	-	-	-	-	-	0.0000	-	0.00632		J	
7 ro04004a	cal 2.50	AC-5-ASA	2.76	33704	2024	6.95	1954948	82708	0.0245	1.0	2.49	-0.5	J	
8	cal 1.00	5-ASA	4.95	18973	1017				0.0123		1.00	0.4	J	
9 ro04005a	cal 6.25	AC-5-ASA	2.77	79998	4782	6.96	1845852	78083	0.0612	1.0	6.23	-0.4	J	
10	cal 2.50	5-ASA	4.96	47587	2426				0.0311		2.53	1.0	J	
11 ro04006a	cal 12.5	AC-5-ASA	2.77	184540	10203	6.93	1959039	83035	0.1229	1.0	12.5	-0.1	J	
12	cal 5.00	5-ASA	4.96	98803	5086				0.0613		4.97	-0.5	J	
13 ro04007a	cal 25.0	AC-5-ASA	2.76	367361	21137	6.95	1993394	84423	0.2504	1.0	25.5	1.8	J	
14	cal 10.0	5-ASA	4.96	202287	10349				0.1226		9.95	-0.5	J	
15 ro04008a	cal 62.5	AC-5-ASA	2.75	872066	51750	6.89	1995195	84562	0.6120	1.0	62.2	-0.4	J	
16	cal 25.0	5-ASA	4.94	494858	25753				0.3045		24.7	-1.2	J	
17 ro04009a	cal 125	AC-5-ASA	2.75	1787784	104513	6.97	2032549	85689	1.2197	1.0	124	-0.8	J	
18	cal 50.0	5-ASA	4.96	1045971	52640				0.6143		49.8	-0.3	J	
19 ro04010a	cal 250	AC-5-ASA	2.75	2751271	162216	6.94	1535522	65686	2.4636	1.0	251	0.4	J	
20	cal 100	5-ASA	4.96	1631033	81657				1.2431		101	0.8	J	
21 ro04011a	cal 625	AC-5-ASA	2.75	7625330	448198	6.94	1719859	72917	6.1467	1.0	625	0.0	J	
22	cal 250	5-ASA	4.95	4571734	226421				3.1052		252	0.7	J	
23 ro04012a	cal 1250	AC-5-ASA	2.77	8879399	520831	6.95	1017799	42367	12.2933	1.0	1250	0.0	J	
24	cal 500	5-ASA	4.95	5279346	260029				6.1375		498	-0.4	J	
25 ro04013a	w.b.	AC-5-ASA	-	-	-	-	-	-	-	1.0	-		J	
26		5-ASA	-	-	-	-	-	-	-	-	-		J	
27 ro04014a	# 1 / 1 / 0.00h	AC-5-ASA	-	-	-	6.95	2253300	95300	0.0000	1.0	-0.000523		0	
28		5-ASA	-	-	-	-	-	-	0.0000	-	0.00632		0	
29 ro04015a	# 1 / 1 / 1.00h	AC-5-ASA	2.75	697857	41275	6.95	1939366	82346	0.5012	1.0	51.0		J	

The list can be reduced to one analyte by double-clicking in the 'analyte' column if more than 1 analyte were measured in one chromatogram. After second double-click in the 'analyte' column the whole batch list is shown again.

file/sample	name	analyte	rt	area	ht	rt(is)	area(is)	ht(is)	p.h.r.	dl.f	ng/mL	dev [%]	flag	comment
1	ro04001a	sss	AC-5-ASA	2.77	3458351	202247	6.98	1951732	82112	2.4631	1.0	250	J	
3	ro04002a	w.b.	AC-5-ASA	-	-	-	-	-	-	1.0	-	-	J	
5	ro04003a	Blank	AC-5-ASA	-	-	6.97	1991341	84158	0.0000	1.0	-0.000523	J		
7	ro04004a	cal 2.50	AC-5-ASA	2.76	33704	2024	6.95	1954948	82708	0.0245	1.0	2.49	-0.5	J
9	ro04005a	cal 6.25	AC-5-ASA	2.77	79598	4782	6.96	1845852	78083	0.0612	1.0	6.23	-0.4	J
11	ro04006a	cal 12.5	AC-5-ASA	2.77	184540	10203	6.93	1953039	83035	0.1229	1.0	12.5	-0.1	J
13	ro04007a	cal 25.0	AC-5-ASA	2.76	367361	21137	6.95	1993394	84423	0.2504	1.0	25.5	1.8	J
15	ro04008a	cal 62.5	AC-5-ASA	2.75	872066	51760	6.89	1995195	84562	0.6120	1.0	62.2	-0.4	J
17	ro04009a	cal 125	AC-5-ASA	2.75	1787784	104513	6.97	2032549	85689	1.2197	1.0	124	-0.8	J
19	ro04010a	cal 250	AC-5-ASA	2.75	2751271	162216	6.94	1535522	65686	2.4696	1.0	251	0.4	J
21	ro04011a	cal 625	AC-5-ASA	2.75	7625330	448198	6.94	1719859	72917	6.1467	1.0	625	0.0	J
23	ro04012a	cal 1250	AC-5-ASA	2.77	8879399	520831	6.95	1017799	42367	12.2933	1.0	1250	0.0	J
25	ro04013a	w.b.	AC-5-ASA	-	-	-	-	-	-	1.0	-	-	J	
27	ro04014a	# 1 / 1 / 0.00h	AC-5-ASA	-	-	6.95	2253300	95300	0.0000	1.0	-0.000523	0	J	
29	ro04015a	# 1 / 1 / 1.00h	AC-5-ASA	2.75	697857	41275	6.95	1939366	82346	0.5012	1.0	51.0	J	
31	ro04016a	# 1 / 1 / 2.00h	AC-5-ASA	2.77	1326328	79377	6.97	1889993	80661	0.9841	1.0	100	J	
33	ro04017a	# 1 / 1 / 3.00h	AC-5-ASA	2.76	2146983	127207	6.97	2250208	95636	1.3315	1.0	135	J	
35	ro04018a	# 1 / 1 / 4.00h	AC-5-ASA	2.77	2259397	134718	6.99	2319572	98761	1.3641	1.0	139	J	
37	ro04019a	# 1 / 1 / 5.00h	AC-5-ASA	2.76	2653472	158542	6.98	2130515	90555	1.7508	1.0	178	J	
39	ro04020a	# 1 / 1 / 6.00h	AC-5-ASA	2.77	2315585	138637	6.99	2139429	91203	1.5201	1.0	155	J	
41	ro04021a	# 1 / 1 / 7.00h	AC-5-ASA	2.76	2663665	158818	6.98	2031510	86135	1.8438	1.0	187	J	
43	ro04022a	# 1 / 1 / 8.00h	AC-5-ASA	2.77	2893913	173513	7.00	1832325	77926	2.2295	1.0	227	J	
45	ro04023a	# 1 / 1 / 9.00h	AC-5-ASA	2.76	2481021	148016	6.98	2106563	89324	1.6571	1.0	168	J	

### Editing Sample Names

After double-clicking (or ENTER) in the 'name' column a dialog appears for entering the new sample name. Times are given in the format " \_\_d\_\_h\_\_m". Expressions such as 2w, 1d2h30m, 5d10m, 20m, 1d etc. are valid. If a number is not followed by a letter the system automatically assumes a h.

Choosing CHANGE ALL TIMES FROM... changes all time information for all SUB samples to the new value. The example below shows that all SUB samples for the current analyte in the project which up to now have the time of 3h will be assigned the new time of 1d12h.

Line No.:	Smpl Type	Subject	Period	Time	Dil.Factor
33	SUB	1	1	3	1.0
	SUB	1	1	1d12h	1

change all times from 3 to 1d12h!

Cancel Ok

The key combination CTRL-first letter of sample type offers a quick way to change the sample type. CTRL-C changes sample type to CAL, CTRL-Q to QC, CTRL-V to VAL, CTRL-S to SUB etc.

## Editing Chromatogram Flags

Die Chromatogram flags indicate acceptance of the results for a specific analyte in a chromatogram. They can be directly entered into the column via the keyboard (only entries valid for the current sample type will be accepted), or the user can 'page' through the possible flags by double-clicking, or the flags are set via the context menu (right mouse button). If several analytes were measured in a chromatogram the flags can be changed for all analytes of the respective chromatogram at the same time by keeping the CTRL key pressed while changing the flag. If the new or old flag is 'D', it will be set for all analytes of the respective chromatogram anyway (of course!).

### Analysts-Flags

Flag	CAL/QCS	VAL	SUB	Explanation
Y	+	+	+	everything OK
N	+	+	+	chromatography error (e.g. peak interferences...), instrument error <b>must be re-assayed</b> (for SUB)
S	+	+	+	chromatographically Ok internal standard not accepted (peak too small/too large) <b>must be re-assayed</b> (for SUB)
X	+	+	+	chromatogram will be ignored in results
D	+	+	+	sample destroyed during work-up  <b>must be re-assayed (only SUB)</b>
E	+	+		chromatographically Ok will not be used for statistic as outside the acceptance limits
R			+	chromatographically Ok, but value not plausible, therefore, marked to re-assay <b>should be re-assayed</b>
C			+	chromatographically Ok, but selected by client for re-assay <b>must be re-assayed</b>
V			+	A "dummy" SUB sample is marked with V (Volume) to report that sample was received but assay/reassay is not possible due to lack of sample material (empty tube) <b>should be re-assayed if in the next sample shipment</b>
I			+	A "dummy" SUB sample is marked with I (Ignore) to report the sample as "Sample received, but not analyzed (NOA)" <b>must not be re-assayed</b>

DBLABCAL sets its own chromatogram flags for chromatograms from unknown samples (SUB) that were accepted by the analyst (with flag Y), and only for these:

**Databases-Flags**

Flag	Explanation
0	chromatographically Ok, calculated concentration <LOQ
<	calculated concentration less than the lowest CAL value in the current batch, but higher than the project LOQ (can happen when the lowest CAL value is not used)
>	calculated concentration higher than the highest CAL values in the current batch
+	calculated concentration higher than the LOQ found in a pre-dose sample (SUB sample with time<=0)
#	double assays acceptance criterion (CV value of the double assays) not fulfilled (only at regression model „double assays)

All chromatograms with flags N, S, D, R, C, <, >, + and # are automatically listed under SAMPLES TO BE RE-ANALYZED.

**The analyst must not set a flag to (for instance) ‘N’ or ‘R’ just because the concentration in a pre-dose sample was <LOQ, >ULQ or >LOQ...**


The respective table’s cell back color is gray to indicate that sample name or chromatogram flag were changed by the user (at least once).

It is possible, with SHOW HISTORY (right mouse button) to view the Audit Trail for the chosen changes.


The screenshot shows the dbLabCal V3 software interface. At the top, there are menu options: Study, Results, Batch, Batch Data, Chromatography, db, Language, Audit Info, Extra. Below this, there are fields for 'Study' (HA107), 'Analyte' (Analyte1 / plasma), and 'Model' (y=a+bx). A 'Sequence' section shows '4' and '9 | Analyst01 | history'. A 'Batch List' table is displayed with columns: file/sampleID, name, analyte, rt, area, ht, rt(is), area(is), ht(is), p.h.r., dLI, ng/mL, dev [%], flag, comment. Row 51 is highlighted in red and has a context menu open over it. The context menu includes options: Y OK, N chrm.error, S irregularity of I.S. area/height, E Cal/Concluded, R to be reassayed (not reliable), C to be reassayed (clients request), D sample destroyed, V insuff. sample volume, X chrm. excluded, and show history!

file/sampleID	name	analyte	rt	area	ht	rt(is)	area(is)	ht(is)	p.h.r.	dLI	ng/mL	dev [%]	flag	comment
37	ro04013a	Analyte1	2.76	2653472	150542	6.98	2130515	90555	1.7509	1.0	178		Y	
38		Analyte2	4.98	1099823	53918				0.5994		48.3		Y	
39	ro04020a	Analyte1	2.77	2315895	13637	6.99	2139429	91203	1.5201	1.0	155		Y	
40		Analyte2	4.98	835998	40924				0.4487		36.4		Y	
41	ro04021a	Analyte1	2.76	2663685	158818	6.98	2031510	86135	1.8438	1.0	187		Y	
42		Analyte2	4.98	1204770	59471				0.6904		56.0		Y	
43	ro04022a	Analyte1	2.77	2899913	173513	7.00	1832325	77826	2.2295	1.0	227		Y	
44		Analyte2	4.99	1249964	61757				0.7935		64.4		Y	
45	ro04023a	Analyte1	2.76	2491021	148016	6.98	2106963	89324	1.6571	1.0	168		Y	
46		Analyte2	4.98	917013	40020				0.4460		36.3		Y	
47	ro04024a	Analyte1	2.76	1276978	76261	6.96	1869157	79093	0.9842	1.0	98.0		Y	
48		Analyte2	4.96	207652	10159				0.1284		10.4		Y	
49	ro04025a	Analyte1	2.75	807667	48607	6.95	1882888	79842	0.6088	1.0	61.9		Y	
50		Analyte2	4.96	107580	4997				0.0626		5.08		Y	
51	ro04026a	Analyte1	2.75	595907	36941	6.93	2308960	94721	0.3784	1.0	38.5		Y	
52		Analyte2	4.95	65641	2897				0.0306		2.49		Y	
53	ro04027a	Analyte1	2.75	263147	16630	6.92	2183090	93659	0.1197	1.0	18.3		N	chrm.error
54		Analyte2	4.95	45284	2012				0.0215		1.75		S	irregularity of I.S. area/height
55	ro04028a	Analyte1	2.75	161664	9518	6.93	2149431	92052	0.1034	1.0	10.5		E	Cal/Concluded
56		Analyte2	4.98	21103	999				0.0109		0.887		R	to be reassayed (not reliable)
57	ro04029a	Analyte1	-	-	-	6.92	1333998	56991	0.0000	1.0	-		C	to be reassayed (clients request)
58		Analyte2	2.74	24218	1284				0.0225		2.29		D	sample destroyed
59	ro04030a	Analyte1	-	-	-	-	-	-	-	1.0	-		V	insuff. sample volume
60		Analyte2	-	-	-	-	-	-	-	-	-		X	chrm. excluded
61	ro04031a	qcs 5.00	2.76	77449	4405	6.91	2110645	91026	0.0484	1.0	4.92	-1.6		
62		qcs 2.00	4.95	41937	2266				0.0243		2.03	1.3		
63	ro04032a	Analyte1	2.75	709538	42337	6.93	2034618	87291	0.4850	1.0	49.3	-1.4		

## Chromatography Menu



The user can view retention times, peak areas, peak heights, peak widths and number of theoretical plates for the current analyte and the internal standard (if used). These evaluations primarily allow for a quick review and check whether the integration software has correctly assigned peaks (retention time, relative retention time) and correctly positioned the baseline (peak width, theoretical plates).



Moving the mouse across the graph automatically displays sample information under the current mouse cursor position in the following order: retention time, peak area, peak height, and peak width.

User may „jump“ to the underlying chromatogram with CTRL-double click.

## Queries...

This offers the chance to select only specific samples types and/or concentrations for the chromatographic evaluation; to view data for all peaks of the chromatogram (instead of only the current analyte) etc...

**Chromatography/Queries**

y-Axes:

- Retention Times
- Peak Areas
- Peak Heights**
- Peak Widths
- theor. Plates
- Concentration

Peak(s):

- Analyte1 / plasma**
- Analyte2 / plasma (is)
- Internal Standard

Sample Type(s):

- CAL
- QCS**
- SUB
- SSS
- VAL (ER: Extracts/Room Temperature)
- VAL (EK: Extracts/Refrigerator (5±3°C))
- VAL (EG: Extracts/Freezer (-20±5°C))
- VAL (ET: Extracts/Deep Freezer (-77±5°C))
- VAL (NR: Matrix/Room Temperature)
- VAL (NK: Matrix/Refrigerator (5±3°C))
- VAL (NG: Matrix/Freezer (-20±5°C))
- VAL (NT: Matrix/Deep Freezer (-77±5°C))
- VAL (PR: Validation Samples)
- VAL (PK: Validation Samples)
- VAL (PG: Pools (freeze/thaw))
- VAL (PT: Validation Samples)
- VAL (BR: Validation Samples)
- VAL (BK: Validation Samples)
- VAL (BG: Validation Samples)
- VAL (BT: Validation Samples)
- VAL (AR: Other Matrix)
- VAL (AK: Validation Samples)
- VAL (AG: Validation Samples)
- VAL (AT: Validation Samples)

Nominal Conc.:

- 2.50
- 5.00
- 50.0
- 500**
- 1000
- 1.00
- 2.00
- 20.0
- 200
- 400

x-Axes:

- no.
- Retention Times
- Peak Areas
- Peak Heights
- Peak Widths
- theor. Plates
- Concentration

Conc.>LOQ only

Ratio Analyte/IS


Cancel Ok

For viewing data to several peaks only the running chromatogram number (no.) is allowed for the x-axis.




Choosing  Conc.>LOQ only displays only data from chromatograms with a calculated concentration that is higher than LOQ. Choosing  Ratio Analyte/IS calculates and displays the size of the y-axis relative to that of the internal standard. Peak area changes to p.a.r. and peak height to p.h.r...

Quality Control Chart data presentation is also possible:  
Nominal concentration (instead of MEAN) and %-deviation lines (instead of sd) are displayed if the number for SD / % is set to  $\geq 5$ .



It is often useful to view peak width as a function of peak height or peak area.


*If, for instance, the following figure would show a peak width of 18 for a peak with a height of 300,000, the respective chromatogram should be closely examined. A peak width of 18 would be acceptable for small peaks.*



*The data on which the graphic display is based can also be exported. This offers a way to export these 'special data' from the data base as well.*

It is also possible to display, print and export various and often very interesting correlations (e.g. between peak width and concentration, or peak height and peak area etc.).

If the same parameters are chosen for the x- and y-axis the percentage distribution density of the data can be viewed.



All graphics as well as the data can be copied at any time into windows clipboard. (CTRL-INS or CTRL-C)


## db Menu

db	QC+QA
	Capacity Browser...
	Long Term Stability Planning
	Attached Users...
	New Password...
	SQL Direct...
	New Admin Password...
	Re-Login required after minutes...
	System AuditTrail
	Edit departments...
	Edit users...
	Edit company/site name...
	Edit database system texts ...
	Edit messages/captions/menus etc...
	<input checked="" type="checkbox"/> Allow *.lca import (Analyst) Force *.rdb for *.lca import
	<input checked="" type="checkbox"/> Allow Empower2 direct import
	<input checked="" type="checkbox"/> Allow *.xls import (Xcalibur)
	<input checked="" type="checkbox"/> Allow *.rep import (ICP-MS)
	<input checked="" type="checkbox"/> Allow *.asc import (Magellan)
	<input checked="" type="checkbox"/> Allow *.csv import (access2)
	<input checked="" type="checkbox"/> Allow *.xls import (SearchLight)
	<input checked="" type="checkbox"/> Allow *.xls import (LabX)
	<input checked="" type="checkbox"/> Allow *.xls import (Chromera)
	Allow *.txt import (SoftMax Pro)
	<input checked="" type="checkbox"/> Allow export to BI ASCII file
	<input checked="" type="checkbox"/> Allow export to HoLaRo ASCII file
	<input checked="" type="checkbox"/> Allow HoLaRo reassay rule
	Allow External access
	<input checked="" type="checkbox"/> Allow Not analyzed flag (NOA)
	<input checked="" type="checkbox"/> Allow Range [%] display
	Allow 2 Languages
	Lock QC/QA items editing
	Allow regression model $y=a+bx+cx^2$
	<input checked="" type="checkbox"/> Allow no calc. model ( $y = x$ )
	<input checked="" type="checkbox"/> Allow regression model $y=a+b.log(x)$
	Allow $1/y$ and $1/y^2$ weghtings
	Max. length of study code...
	Max. sequences...
	Max. lines for results...
	Max. nominal conc. for validation/stability samples...
	Max. chrmdata entries/project (results auto load)...

Most entries (starting with the first line) are only accessible to the dbLabCal administrator.

## Attached users...

This menu item is displayed for department managers and higher authorities. Through this menu the department manager can log off users (by 'clicking away' the name) who, because of incorrectly exiting the program, are still logged into the database and may therefore block the editing of project data by other users.




## SQL Direct...

This menu allows for extracting additional information from the database. Examples: status of all studies in the database, number of samples measured within a specific period of time, with a specific analytical method etc. As the paragraph title indicates this process involves ad-hoc SQL queries of the database. Even though only SELECT queries are allowed they still can keep the database busy for a while, especially if the sql statement was incorrect.

In order to be able create own SQL queries the authorized users must be familiar with both SQL and the structure of the DBLABCAL database.

Some SQL queries are already predefined after installation.



Up to four SQL Direct dialogs can be opened simultaneously.

## Capacity Browser

Der berechtigte Benutzer hat hier einen Überblick über die Ausnutzung der einzelnen Messplätze (Kapazitätsauslastung) und über die Anzahl der gemessenen Proben pro Tag oder pro Monat

## Long Term Stability Planning

This dialog is useful to plan long term stability experiments.

It shows also all stability data results. All results are maintained automatically by dbLabCal. It means, dbLabCal adds the result into this list immediately after respective VAL samples were imported or renamed.

Code	Project	Analyte	Matrix	Project Type	Analyte Stat.	Months Plan	Due Date	Reminder	ValD Done	LTS Done	Condition	Result	Comment	finished
1	BA294	MFA in Serum	MFA	Serum	Validation	data releaser	1	15.07.2002	02.10.2002		MATRIX / GS	5.2		NO
2	BA598	Langzeitstab	15-Keto-PGE0	Plasma	Validation	started	1	11.06.2002	08.06.1999		MATRIX / TK	6.5		NO
3		Langzeitstab	PGE0	Plasma	Validation	started	1	11.06.2002	08.06.1999		MATRIX / TK	6.5		NO
4		Langzeitstab	PGE1	Plasma	Validation	started	1	11.06.2002	08.06.1999		MATRIX / TK	6.5		NO
5	EA004	Amlodipine in	Amlodipine	Plasma	Validation	data releaser	1		08.01.2002		MATRIX / GS	73.0		NO
6	EX011	Amoxicillin in	Amoxicillin	Plasma	Validation	data releaser	1	05.03.2001	05.03.2007		MATRIX / TK	2.8		NO
7	PX037	Clavulanic Ac	Clavulanic Acid	Plasma	Validation	data releaser	1	06.03.2001	06.03.2007		MATRIX / TK	2.8		NO
8	PX064	Interday	Clonidine1997LTS	Plasma	Validation	data releaser	1		12.06.1997		MATRIX / GS	3.0		NO
9	GA329	4-Oxo-/Iso-/	4-oxo-Isotretinoin	Plasma	Routine	data releaser	1	01.12.1998	09.02.1999		MATRIX / GS	2.4		NO
10		4-Oxo-/Iso-/	Isotretinoin	Plasma	Routine	data releaser	1	01.12.1998	09.02.1999		MATRIX / GS	2.4		NO
11		4-Oxo-/Iso-/	Tretinoin	Plasma	Routine	data releaser	1	01.12.1998	09.02.1999		MATRIX / GS	2.4		NO
12	GA497	Trimethoprim	Sulfamethoxazol	Plasma	Routine	data releaser	1	02.09.1998	15.10.1999		MATRIX / TK	13.6		NO
13		Trimethoprim	Trimethoprim	Plasma	Routine	data releaser	1	02.09.1998	15.10.1999		MATRIX / TK	13.6		NO
14		Trimethoprim	Sulfamethoxazol	Urin	Routine	data releaser	1	15.09.1998	15.11.1999		MATRIX / GS	14.5		NO
15		Trimethoprim	Trimethoprim	Urin	Routine	data releaser	1	09.09.1998	25.10.1999		MATRIX / GS	13.7		NO
16	GA526	Methimazole	Methimazole	Plasma	Validation	data releaser	1	22.05.1998	26.09.1998		MATRIX / TK	4.3		NO
17	GA638	Revalidation	Ro 1-8551_Dog	dog plasma	Validation	data releaser	1	12.03.1998	07.09.1998		MATRIX / GS	5.4		NO
18		Revalidation	Ro 1-8551_Rat	rat plasma	Validation	data releaser	1	12.03.1998	07.09.1998		MATRIX / GS	5.4		NO
19		Revalidation	Ro 1-8551_Rat	rat plasma	Validation	data releaser	1	06.03.1998	08.09.1998		MATRIX / GS	6.1		NO
20		Revalidation	Ro 1-8551_Rat	rat plasma	Validation	data releaser	1	06.03.1998	08.09.1998		MATRIX / GS	6.1		NO
21	GA671	Metabolite M	M19	Plasma	Validation	data releaser	1	22.04.1998	22.09.1998		MATRIX / GS	5.3		NO
22	GA702	VALSARTAN	Valsartan_V	Plasma	Validation	data releaser	1	13.10.1998	12.01.1999		MATRIX / GS	3.3		NO
23	GX004	Buspirone in	1-PP	Plasma	Validation	started	1		13.03.2002		MATRIX / GS	39.6		NO
24		Buspirone in	Buspirone	Plasma	Validation	started	1		13.03.2002		MATRIX / GS	39.6		NO
25	GX006	LZST CE frei	Delta9-9DHE	Plasma	Validation	data releaser	1		18.10.2002		MATRIX / GS	6.3		NO
26		LZST CE frei	Equilin	Plasma	Validation	data releaser	1		18.10.2002		MATRIX / GS	6.3		NO
27		LZST CE frei	Estroren	Plasma	Validation	data releaser	1		18.10.2002		MATRIX / GS	6.3		NO
28	GX010	Loratadine/C	Descarboethoxyflora	Plasma	Validation	data releaser	1		04.01.2000		EXTR / RT	1.1		NO
29		Loratadine/C	Loratadine	Plasma	Validation	data releaser	1		04.01.2000		EXTR / RT	1.1		NO
30	GX020	Testosteron in	Testosteron	Plasma	Validation	data releaser	1	13.05.2003	19.05.2003		EXTR / RT	10.3		NO
31	GX039	Gemfibrozil in	Gemfibrozil	Plasma	Validation	data releaser	1	11.01.2007	11.01.2007		MATRIX / GS	3.1		NO
32	GX047	Validierung C	CPA	Plasma	Validation	data releaser	1	29.05.1998	29.05.1998		MATRIX / GS	8.2		NO
33	GX048	Estronsulfate	Estron	Plasma	Validation	data releaser	1	10.03.1998	10.03.1998		MATRIX / GS	1.1		NO
34	HA058	Peindopril/P	Peindopril	Plasma	Validation	data releaser	1	27.04.1998	17.06.1998		MATRIX / GS	1.7		NO
35		Peindopril/P	Peindoprilat	Plasma	Validation	data releaser	1	27.04.1998	17.06.1998		MATRIX / GS	1.7		NO

## Language

Here, the user can toggle at any time between the two pre-defined languages. In such a switch the data of the current project are reloaded and redisplayed in the newly selected language.

## New Password...

The user can change his/her Oracle password at any time through this menu.

*The following menu options are only accessible to the DBLABCAL administrator*

**New Admin Password...**

The administrator can change the Oracle admin password at any time through this menu.

**Re-Login required after minutes...**

If the user logged into DBLABCAL is different from that logged into the PC/network and if over the period of time specified in this option no activity is registered, the user logged into DBLABCAL will automatically be logged off for security reasons.

**Edit company/site name...**

The name of the company and/or site as it should appear on the DBLABCAL printouts is entered through this option.

Edit departments...

The screenshot displays a 'Departments' dialog box with a blue title bar. At the top, it shows 'ID: 0', 'Name: Bioanalytics', and a checked 'has instruments' checkbox. There are 'New' and 'Exit' buttons to the right. The main area contains a table with two columns of data, each with an 'in use' checkbox, an 'Instrument' column, and a 'Comment' column.

in use	Instrument	Comment	in use	Instrument	Comment
<input checked="" type="checkbox"/>	LCMS-01	API 5500+Analyst 1.5.1	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	LCMS-02	API 5000+Analyst 1.4.2	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	LCMS-03	API 4000+Analyst 1.4.2	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	GCMS-01	DSQ+Xcalibur 2.0	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	GCMS-02	DSQII+Xcalibur 2.0	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	HPLC-01	HPLC-FL+Empower	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	HPLC-02	HPLC-EC+Empower	<input type="checkbox"/>		
<input checked="" type="checkbox"/>	ELISA	ELISA+Magellan	<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		
<input type="checkbox"/>			<input type="checkbox"/>		

At the bottom of the dialog are 'Cancel' and 'Save' buttons on the right, and 'First', 'Previous', 'Next', and 'Last' buttons on the left.

The department names are entered through this dialog. If the respective department belongs to QM, for instance, and therefore it is not an analytical department, the check box HAS INSTRUMENTS will be deactivated.



**Edit users...**

The screenshot shows a dialog box titled "Users" with the following fields and controls:

- ID: 254 (with "New" and "Edit" buttons)
- User ID: VALANALYST (with an "account" checkbox)
- Name: VAL-Analyst (with a "locked" checkbox)
- Password for Oracle Login: (empty text field)
- Department: VALDEPT (dropdown menu)
- Authorization: Analyst (dropdown menu) with an "active" checkbox
- Checkboxes:  'W' flag allowed,  QA-Reviewer,  external user,  released data only
- SQL Direct: Select Query Only (dropdown menu)
- Capacity Browser: Own Department (dropdown menu)
- Buttons: Find!, Cancel, Save!, First, Previous, Next, Last

ACTIVE is activated if the user is an analyst. This ensures that this analyst is displayed in the lists of analysts during the editing of a project or the import of batches. DBLABCAL does not delete any users (and also no departments). If a user does not work with DBLABCAL, only ACTIVE will be reset and/or the Oracle user locked in Oracle database.

UserID is the user name with which the user logs into the network. The name appears in dbLabCal in the lists of analysts and also in the status line. Permission levels are: ReadOnly, Analyst, PI/Study Director, Dept. Manager and Administrator. SQL Direct is used for specifying whether a specific user sees the SQL Direct menu at all and whether he/she can write SQL queries. "R' flag allowed" is always enabled for users with access authority of PI/study director and higher.

„QA-Reviewer“ is a ReadOnly user with the additional permission to edit the QC/QA dialog.

Check boxes ACCOUNT and LOCKED show current Oracle status of the user. Does user exist in Oracle DB? Is user locked by Oracle DB?

**Edit messages/captions/menus etc...**

**Messages**

ID: 10           

Description: Abbrechen

Cancel

---

Tab : >>    New line : //    Parameters : \_1\_, \_2\_, \_3\_ or \_4\_

german: Abbrechen

english: Cancel

In this dialog all program texts, annotations etc. can be changed/edited. Any changes are immediately visible to the administrator, all other users will see them only after DBLABCAL is restarted.

**Edit database system texts...**

**System Text**

ID: 1           

Description: Calibration sample

Cancel

---

Tab : >>    New line : //    Parameters : \_1\_, \_2\_, \_3\_ or \_4\_

german: CAL

english: CAL

This dialog allows for defining the sample type (CAL, QCS, SUB, etc.), codes for matrix and temperature (EXTR, RT, GS, etc.), codes for chromatogram flags, codes for batch status, project type, batch type, readings, audit/trail texts etc. Any changes are shown after restarting DBLABCAL.

### **Allow...**

This dialog allows for specifying which menu options are visible and thereby defines the overall user functionality of DBLABCAL.

### **Max...**

If desired, this item allows for extending the database limits.

### **Menü Extra**

This menu is visible if set-up by the administrator. It allows to start any other program including extensions of DBLABCAL. Purpose and use of such extra programs must be described separately.

---

## **Tips and Tricks**

The change in the mode of working with analytical data when going from a 'paper' office to an 'Excel' office using the database is reiterated. The major points are as follows:

### **Naming of samples:**

This must be well thought out since the results are 'tightly connected' to the sample names. It is possible to change the samples names later in the database (batch status, audit trail), but this is an elaborate, time-consuming process involving each time a complete recalculation of batch data and update of the project results.

### **Setting of chromatogram flags, batch status and project status**

Chromatogram flags, batch status and project status are at least as important as the analytical values and have well defined meanings and conbatches.

*If sample names and/or chromatogram flags are changed in the database later on, the user must not forget to make the same changes on the chromatograms as well.*

The data received from the database is still as correct or as incorrect as the data written into the database. The easier the 'not quite correct' data are passed on the more important it is that the raw data imported into the database are 100% checked.



The data path from unit to report is as follows:

- 1) Generation of chromatograms (**raw data**) at the unit  
ANALYST, EMPOWER2...
- 2) Summarizing the chromatographic results of a batch in an ASCII file  
EXPORT-OPTION DER CHROMATOGRAPHIE-SOFTWARE, QUANCAL, MANCAL
- 3) If required (when units are not connected to network) transferring the ASCII file generated by copying the file onto data server
- 4) **Importing** the ASCII file into the database  
From here on ASCII files are not needed. Changes to or deletion of the files have no impact whatsoever on the data in the database.
- 5) Exporting the data from the database into the report  
DBLABCAL, WORD, EXCEL

## The number are „unrealistic high“

dbLabCal works at best with windows international setting (numbers,currency) set to US. It can have problems with the representation of numbers if both decimal separators (period and comma) are selected in the number/currency format.

Solution: check settings in workplace – control panel – country settings/currency format. The separator for groups of numbers should be set to blank.



### **A chromatogram imported into the database must be reintegrated**

Since the data are imported batch-wise into the data base any modification in the data of a chromatogram necessitates a full "import" process, that is:

If the user notices in the database (?!) that the integration software generated incorrect peak assignments or baselines, the following actions are required:

- 1) (Manually) re-integrate the chromatogram with the integration software
- 2) Ensure that the new results can be further processed with QUANCAL, MANCAL (i.e. perform "save as ASCII" if required)
- 3) Import the new file into the database  
If the original batch still shows the 'not set yet' status, this is no problem;  
if the original batch shows the 'batch accepted' status it must be set to 'excluded' beforehand since an 'accepted' batch must not be overwritten

### **The program incorrectly generates the message: "locked by XYZ"**

XYZ could not log off the database because a malfunction terminated the program. XYZ should load the respective project and re-exit it by selecting another project or terminating the program.

The PI/study director can also manually log off the user (DB | ATTACHED USERS...).

### **Changing readings or weightings in a routine project**

Readings or weightings can only be changed in a validation project. Therefore the project type must first be changed to validation. Then a new selection can be made via the menus PROJECT | WEIGHTING or PROJECT | READINGS. All data for all analytes that were measured in a chromatogram will be recalculated based on the new settings. Don't forget afterwards to reset the project type back to routine.

**Subject samples must be measured twice (e.g. measuring protein binding)**

In order that the program does not output the twice-measured samples as re-assays they must not carry the same name. This is best achieved by changing the period, e.g.:

subject 1	period 1	4.00 h	plasma
subject 1	period 2	4.00 h	plasma
subject 1	period 3	4.00 h	actually period 1, buffer
subject 1	period 4	4.00 h	actually period 2, buffer

**User wants to print batches, but not nothing prints**

This is possible with batches still having the ,not set yet' status and the user forgot to deselect the option ONLY 'OK' BATCHES.

**Printing figures with subject concentration curves**

DBLABCAL does not make provisions for printing subject concentration curves. However, the figure can be exported into another program and printed from within there.

- 1) Press CTRL-INSERT (or CTRL-C), select GRAPHIC, or right mouse over the picture
- 2) Change into the target program
- 3) In the Edit menu of the target program select PASTE or a similar expression, or press SHIFT-INSERT (or CTRL-V).

You also can generate your own presentation of the concentration curve by exporting the (selected) data from DBLABCAL into, for instance, Excel.

### **Saving data in an ASCII file**

An ASCII file with TABs is best be generated with Windows Notepad as follows:

- 1) Export data into Notepad - either with CTRL-INSERT (or CTRL-C),  
or right mouse button at marked data
- 2) Select PASTE in the Edit menu of Notepad,  
or press SHIFT-INSERT (or CTRL-V)
- 3) Save file

### **In subjects, some time values shown twice, with 'missing' or with value...**

This happens if some (one is already enough!) samples were assigned only slightly different times. Example: A sample was named SUBxxDGxxTP 0.17h, and later perhaps SUBxxDGxxTP 10m. Since 0.17h is basically 10.2min DBLABCAL differentiates between the two times even though (rounded) each time 10m (or 0.17h) is shown. Change in this case for one SUB sample with time 0.17h the time to 10m and select ALL TIMES FROM... in the dialog SAMPLE NAME (see page 67)

**After selecting the validation samples results the following message displays:  
"No validation 0 values found..., check sample names and/or dates..."**

This exactly means what it says. In most cases the database lacks information about the date "start of batch". Best course of action is to enter the missing date in RESULTS | DATE OF ANALYSES. If the database does contain all dates to the start of batch the reason may be found in missing information about matrix and/or temperature in the sample names for the validation 0 results. You should load the corresponding batch(s) and check for correct sample names.



---

**In a project (the same) analytes are defined in different matrices. Samples with different matrices were simultaneously measured in a batch. However, DBLABCAL only allows a sequential import into one specific analyte/matrix pair defined for a project**

The file to be imported contains the information about the analyte(s) measured in the batch and the (one) matrix. During the import, DBLABCAL checks this information and compares it with that in the current project. If DBLABCAL finds a match, the batch data will be imported into the database.

In such case, the ASCII file has to be imported twice:

then imports the file and manually assigns the first analyte/matrix pair, and then re-imports the same file and manually assigns into the next analyte/matrix pair, etc. The chromatogram flags of the chromatograms with the 'other' matrix will then be set in the database to 'X' – or will be renamed (DIV).

*Describing this situation is much more complicated than actually doing it. Therefore an example:*

The following are defined in the project: - analyte1 in rat plasma

- analyte1 in dog plasma

The following were measured in the batch: - calibration curve (analyte1 in whatever)

- 50 samples of analyte1 in rat plasma

- 10 samples of analyte1 in dog plasma

In generating the batch file (see the respective program) analyte1 will be defined in matrix 'XXX'

- During import DBLABCAL does not recognize analyte1 in matrix 'XXX' and prompts for a manual peak assignment.
- The user selects 'analyte1 in rat plasma'
- After the batch has been imported the user sets the chromatogram flags of the 10 chromatograms of analyte1 in dog plasma to 'X' – or renames them into DIV.
- The batch file is re-imported. DBLABCAL again prompts for a manual peak assignment.
- This time the user selects 'analyte1 in dog plasma'
- The chromatogram flags of the 50 chromatograms of analyte1 in rat plasma are now set to 'X' - or they are renamed into DIV.
- Finished

### **Program limits**

The program limits are chosen such that the vast majority of studies can be processed without undue demands on the PC's resources.

The program limits are set in a table of the database and can only be modified by the administrator. If during the processing of a project the program limits are exceeded, an appropriate message will display.

---

## **Errors and Error Messages**

Apart from the case mentioned above and a number of other questions, notes and security queries, there are three different kinds of error messages:


- 1) The message clearly states what happened where.  
Example: "Error when drawing the dots...", or "Problems with OLE". Action aborted...". This type of error is not really critical even though "something" is wrong with the data.
- 2) The message starts with "Unexpected error...".  
This is more critical. Even though the program is still stable it is not sure for how long. All "Unexpected errors" are logged in the database so that the administrator can retrieve all required information. The administrator should be informed about any occurrences of such an error.
- 3) Short error message.  
This means that right after clicking on OK the program is gone! Probably due to serious problems with computer HW or SW. Re-boot the PC and try again.  
In case the error message can be reproduced make sure to note all of the following information and pass it on to the programmer:  
Project, analyte, batch, and last action performed.  
Since it is highly unlikely that the reason for such a situation is a bug :- ) the user should always first check whether the server may have been unavailable for a short period of time or whether the computer may be defect...

**Appendix 1: HoLaRo- and. BI- ASCII-Files**


This menu option is only visible if enabled by the administrator and if the project has been released.

HoLaRo ASCII is a well-defined ASCII format from Hoffman-LaRoche Company for transferring analytical results. Similarly, Boehringer-Ingelheim Co. has defined an own ASCII format. Data with this ASCII format can be generated via the following dialogs.


Project load, Project Menu -> Results for BI...



HoLaRo-Dialog:




BI-Dialog:



IMPORT DB allows for first loading the data contained in DBLABCAL into the table. It is also possible to load data from an ASCII file (text with TAB separators) with IMPORT ASCII. The data can be edited in the above dialog (or in Excel). The data can be exchanged at any time via the clipboard (or an ASCII file) between Excel and DBLABCAL for editing. Name and path of the file(s) can be entered directly or via '...'.

dbLabCal loads „as good as possible“ the data from the database into the dialog



„as good as possible“ means:

BI:

dbLabCal	BI-ASCII-File-Column
Subject	PTNO
Period	ACTEVENT
Time	PTM und PTMSP

HoLaRo:

dbLabCal	HoLaRo-ASCII-File-Column
subPATIENT_ID_SAMPLE_ID_SAMPLENO	PATIENT_ID
subPATIENT_ID	
subPATIENT_IDdg1tpSAMPLENO	
Period: N/A	SAMPLE_ID
Time : N/A	SAMPLENO

In this dialog all data could be modified (in this dialog ONLY, not in the database).

The dialog data need very likely to be edited according to the DTA (data transmission agreement)

i.e.. „TIOTROPIUM“ to „BA 679“ or „LC-MS/MS“ to „HPLC-MS/MS“ etc...

That's simple and fast with following keys or key combinations:

#### Edit a CELL

Double click or ENTER in cell	->	opens cell for editing
ENTER in opened cell	->	closes cell with new content
ESCAPE in opened cell	->	closes cell with original content

#### SEARCH&REPLACE

CTRL-D in opened cell	->	fills cells below current cell with current cell content
SHIFT-CTRL-D in opened cell	->	fills EMPTY cells below current cell with current cell content
CTRL-R in opened cell	->	replaces all cells below current cell with same content as the current cells original content with the current (visible) content

#### Edit ROW

DEL / INS (cell closed)	->	row will be deleted / empty row inserted
-------------------------	----	--

CTRL-C / CTRL -INS / right mouse + Copy with selected cells

	->	copy into clipboard
CTRL-V oder SHIFT-INS	->	copy from clipboard

It is possible to exchange the data between this dialogs and e.g. Excel at any time via windows clipboard.

Choose file name according toDTA

Lastly, the file is generated in ASCII format with WRITE TO HOLARO or WRITE TO BI!

**IMPORTANT!** : All exported data in the ASCII file must be checked with DBLABCALs data!

---

**Appendix 2: PE SCIEX LC-MS (Analyst)**

Data are exported from the PE SCIEX Software ANALYST (V 1.1 – 1.5.2 were tested) as follows.

In ANALYST start the QUANTITATION WIZARD, select QUANTITATE.

After selecting the batch file (\*.wiff) and the chromatograms in the wizard the RESULTS TABLE will be displayed.

The following columns **must** be defined in this RESULTS TABLE (if not yet done...).

<b>Sample Name</b>	Analyte Peak Name	Analyte Concentration	...
...	Analyte Retention Time	Analyte Peak Area	Analyte Peak Height
...	IS Retention Time	IS Peak Area	IS Peak Height
...	Acquisition Date	Sample ID	Dilution Factor

The **first column must** be SAMPLE NAME -it will be shown in Analyst by default as first column.

The order of the other columns is arbitrary. Additional columns may be defined in the RESULT TABLE.

**The RESULT TABLE must be sorted according to chromatograms – e.g. in the measurement sequence.**

After selecting EXPORT... in the FILE... menu the system prompts for a name for the ASCII file to be exported. The user should enter a name that contains any two letters and not more than three numbers (=batch number). File name extension should be \*.LCA. This 'tells' DBLABCAL that the file is a ASCII file generated by Analyst. (Example: sq01.lca, we7.lca, ab115.lca...)

If there is no Sample ID found in the Result Table, dbLabCal creates a „file name/Sample ID“ from first 5 characters of the lca file and the vial position.

The \*.lca ASCII file is imported into dbLabCal via menu BATCH IMPORT...

### **Appendix 3: Empower2**

In order to export data, it is required that following custom fields are defined . These user-defined custom fields are defined in Empower2:

SType: type of sample, see table below  
 Info\_Sub: subject number for SUB, concentration for CAL, QC or VAL samples  
 Period: period number  
 STime: subject sampling time or storage duration for validation/stability samples  
 Matrix: matrix group used in validation/stability samples  
 Temp: temperature group for validation/stability samples

#### Empower2 naming convention

<b>SType</b>	<b>Info_Sub</b>	<b>Period</b>	<b>STime</b>	<b>Matrix</b>	<b>Temp</b>
cal	nom. concentration *)				
qc	nom. concentration				
sub	number or name	Period	sampling time in h or _d_h_m		
eqc	number or name				
val	nom. Concentration		Dauer in h oder _d_h_m	E Extracts N Matrix P freeze/thaw A A group B B group	R (RT) K (KS) G (GS) T (TK)
blank	e.g. with or without ISor batch number, etc				
diverse	Name				
sss					
w.b.					

\*) first dbLabCal „tries“ to get nom. concentration for CAL from the Amount value in Empower2. If Amount was not used in Empower2, dbLabCal extracts the nominal concentration from the Info\_Sub field.

Two possibilities for combining the fields Info\_Sub, Period and STime for SUB samples exist.

Either

Info\_Sub and Period are entered as **numbers** and STime in hours or as \_d\_h\_m, then dbLabCal interprets the entries as “subject/period/time“

or

if even one entry is **not a number** (or \_d\_h\_m for STime) dbLabCal interprets the entries as “subject“ with the fields Info\_Sub, Period and STime being connected by an underscore.

Apart from the sample name field and the custom fields the following data must be exported as well: vial, dilution, and (of course) retention times, peak areas and peak heights.

Data from Empower2 are imported directly without using any ASCII file with a special dbLabCal's dialog „Import from Empower2...“.

User has to login into Empower2 entering Empower2 userID and password. Current project is selected in the drop down box if found in Empower2.

User selects SampleSet (or SampleSets), then chromatograms to be imported. Selected data are imported after click on **Ok** button.

**Import from Empower2...**

Database: **NEUEMP2** | UserID/PWD: **vagadaym / \*\*\*\*\*** | **login** | Projects: **TA348\_Acyclovir** | **Launch Empower**

Show Sample Sets for selected Project: 1 sample set(s) selected | Show latest Results for selected Sample Sets: 90 chrn(s) selected

SampleSetName	SampleSetStartDate	SampleSetFinishDate	SystemName	Vial	SType	Info_Sub	Period	STime	Temp	Matr	Dilution	DateAcquired	ResultId
run17_unit03	24-Oct-2009 06:11:33	25-Oct-2009 05:42:15	MP03	1	k-ref						1,0000	24-Oct-2009 06:28:27	7754
run18_unit07	23-Oct-2009 21:25:20	24-Oct-2009 23:07:44	MP07	2	w.b.						1,0000	24-Oct-2009 06:44:49	7674
run16_unit03	23-Oct-2009 05:40:51	24-Oct-2009 06:11:32	MP03	3	blank						1,0000	24-Oct-2009 07:01:11	7675
run15_unit07	22-Oct-2009 19:41:58	23-Oct-2009 21:25:19	MP07	4	blank + IS						1,0000	24-Oct-2009 07:17:31	7676
run14_unit03	22-Oct-2009 07:44:45	23-Oct-2009 01:27:26	MP03	5	cal	10.0 ng/mL					1,0000	24-Oct-2009 07:33:50	7672
run13_unit07	21-Oct-2009 19:43:38	22-Oct-2009 19:41:57	MP07	6	cal	20.0 ng/mL					1,0000	24-Oct-2009 07:50:10	7664
run12_unit03	21-Oct-2009 06:40:23	22-Oct-2009 07:44:44	MP03	7	cal	40.0 ng/mL					1,0000	24-Oct-2009 08:06:32	7665
run11_unit07	20-Oct-2009 18:00:17	21-Oct-2009 19:43:37	MP07	8	cal	80.0 ng/mL					1,0000	24-Oct-2009 08:22:52	7666
run10_unit03	20-Oct-2009 06:24:10	21-Oct-2009 06:40:22	MP03	9	cal	100 ng/mL					1,0000	24-Oct-2009 08:39:12	7667
				10	cal	200 ng/mL					1,0000	24-Oct-2009 08:55:33	7668
				11	cal	400 ng/mL					1,0000	24-Oct-2009 09:11:53	7669
				12	cal	800 ng/mL					1,0000	24-Oct-2009 09:28:13	7670
				13	cal	1000 ng/mL					1,0000	24-Oct-2009 09:44:34	7671
				14	w.b.						1,0000	24-Oct-2009 10:00:54	7677
				15	sub	23	1	0h			1,0000	24-Oct-2009 10:17:14	7755
				16	sub	23	1	0.5h			1,0000	24-Oct-2009 10:33:33	7756
				17	sub	23	1	1.0h			1,0000	24-Oct-2009 10:49:54	7757
				18	sub	23	1	1.5h			1,0000	24-Oct-2009 11:06:14	7681
				19	sub	23	1	2.0h			1,0000	24-Oct-2009 11:22:34	7758
				20	sub	23	1	2.5h			1,0000	24-Oct-2009 11:38:54	7759
				21	sub	23	1	3.0h			1,0000	24-Oct-2009 11:55:15	7760
				22	sub	23	1	3.5h			1,0000	24-Oct-2009 12:11:34	7761
				23	sub	23	1	4.0h			1,0000	24-Oct-2009 12:27:54	7762
				24	sub	23	1	5.0h			1,0000	24-Oct-2009 12:44:14	7763
				25	sub	23	1	6.0h			1,0000	24-Oct-2009 13:00:34	7764
				26	sub	23	1	8.0h			1,0000	24-Oct-2009 13:16:54	7765
				27	sub	23	1	10h			1,0000	24-Oct-2009 13:33:14	7766
				28	sub	23	1	12h			1,0000	24-Oct-2009 13:49:37	7767
				29	sub	23	1	16h			1,0000	24-Oct-2009 14:05:59	7768
				30	sub	23	1	24h			1,0000	24-Oct-2009 14:22:20	7693
				31	w.b.						1,0000	24-Oct-2009 14:38:40	7694
				32	sub	23	2	0h			1,0000	24-Oct-2009 14:55:00	7769
				33	sub	23	2	0.5h			1,0000	24-Oct-2009 15:11:20	7770

Unit: **MP03** | HPLC MP03 | Analyst: **\Wick** | Batch No: **17** |  overwrite

**Confirm to overwrite batch!**  
**status - batch status not set yet**  
 extracted  
 Comment to Sequence: **sample set name=run17\_unit03 (id 7078)**

Project: **TA348\_Acyclovir** | Unit: **MP03** | Batch No: **17** | Comment: **extracted**  
 started: **24-10-2009 06:28:00** | finished: **25-10-2009 05:26:00** | ChrmCount: **89**

Peak Assignment: **Acyclovir / Plasma** | Peak 1: **Acyclovir**

**Cancel** | **Ok**



#### **Appendix 4: Xcalibur**

Start QUANBROWSER in Xcalibur (Version 1.4 was tested). Load the sequence file (\*.slr) into QUANBROWSER from the data folder (e.g.: c:\Xcalibur\data\PROJECT\BATCH\ )

Choose EXPORT DATA TO EXCEL in FILE menu of QuanBrowser, then choose EXPORT LONG EXCEL REPORT.

The file name should be a combination of 2 characters and maximally 3 numbers (for the batch number). File extension is always \*.xls.

The \*.xls file is imported into dbLabCal via menu BATCH IMPORT...

## Appendix 5: SoftMax Pro

Create the export file (\*.txt, plain ASCII with TABs) as described in SoftMax Pro. File structure should be as in the following example:

COL	1	2	3	4	5	6	7	8	9	10	11
ROW											
1	##BLOCKS= 4										
2	Group:	DONT	1								
3	Sample	Wells	Sample#	Values	MeanVa						
4	SSS01	A3	1	0.327	0.328						
5		A4		0.328							
6											
7	~End										
8	Group:	SUB	1								
9	Sample	Wells	Sample#	Values	MeanVa	BackCal	Mean Conc.				
10	SUB01DGITP1	D1	1	0.108	0.112	5.78	5.42				
11		D2		0.117			6.14				
12	SUB01DGITP2	D4	2	0.146	0.146	8.45	8.44				
13		E3		0.146			8.46				
14	SUB01DGITP3	D5	3	0.12	0.121	6.45	6.43				
15		D6		0.121			6.47				
16	SUB01DGITP4	D7	4	0.089	0.089	3.96	3.95				
17		D8		0.09			3.97				
18	SUB01DGITP5	D9	5	0.049	0.049	0.71	0.71				
19		D10		0.049			0.71				
20											
21	~End										
22	Group:	QCS	1								
23	Sample	Wells	Concent	Values	MeanVa	BackCal	% Bias	Mean			
24	QCS 15	C1	15	0.225	0.223	14.81	-1.3	14.7			
25		C2		0.222		14.6	-2.6				
26		C7		0.22		14.43	-3.8				
27		C8		0.226		14.89	-0.7				
28	QCS 45	C3	45	0.588	0.585	44.83	-0.4	44.6			
29		C4		0.587		44.76	-0.5				
30		C9		0.573		43.55	-3.2				
31		C10		0.591		45.11	0.3				
32	QCS 80	C5	80	0.961	0.972	77.69	-2.9	78.6			
33		C6		0.992		80.45	0.6				
34		C11		0.97		78.51	-1.9				
35		C12		0.963		77.84	-2.7				
36											
37	~End										
38	Group:	CAL	1								
39	Sample	Wells	Mean Co	Concent	BackCal	Mean Co	Values	MeanV	Std.Dv	CV%	% bias
40	Cal 5	B1	5.02	5	5.07	5.02	0.103	0.1	0	0.8	1.4
41		B2			4.98		0.102				-0.4
42	Cal 10	B3	9.9	10	9.92	9.9	0.164	0.16	0	0.3	-0.8
43		B4			9.88		0.164				-1.2
44	Cal 25	B5	25.22	25	25.31	25.22	0.354	0.35	0	0.4	1.2
45		B6			25.13		0.352				0.5
46	Cal 50	B7	49.7	50	50.17	49.7	0.65	0.64	0.01	1.2	0.3
47		B8			49.23		0.639				-1.5
48	Cal 75	B9	75.25	75	75.75	75.25	0.94	0.93	0.01	0.8	1
49		B10			74.76		0.929				-0.3
50	cal 100	B11	99.34	100	99.37	99.34	1.195	1.2	0.01	0.7	-0.6
51		B12			100.5		1.207				0.5
52											
53	Smallest standard value:	0.103									
54											
55	Largest standard value:	1.201									
56											
57	~End										

- use dbLabCals naming conventions for the Group name (CAL, QCS, SUB, VAL, SSS, DIV) if you want import data of that group (if you don't want to import data from a group, just name the group with some text unknown to dbLabCal)
- use dbLabCals naming conventions for the sample name (CAL, QCS, SUB, VAL, SSS, DIV) it is sufficient to name the samples only once, if the following row has data of the same sample, as shown above

- define the data columns for each sample type in the softmaxpro.col file in the dbLabCals working folder

softmaxpro.col file – example:

```
[default]
DataColCAL=5
DataColQCS=6
DataColSUB=7
DataColSUBDILF=-1
DataColVAL=6
DataColKON=6
DataColDIV=5
```

- sample name is in the first column for all sample types

The \*.txt file is imported into dbLabCal via menu BATCH IMPORT...

All individual data (not only the mean values!) will be imported.

It means:

- there is NO regression calculation in dbLabCal, and therefore, NO concentration calculation for QCS, SUB, VAL and SSS samples
- the user will have to make choice on final SUB results for ALL samples since for ALL SUB samples at least two individual results were imported.

**Appendix 6: Magellan (ELISA Reader)**

Magellan data (Magellan V 3.11) with extension “asc” are saved automatically in pre-defined folder on the file server.

The \*.asc file is imported into dbLabCal via menu BATCH IMPORT...

**Important: asc file import is possible only if current project uses regression model NO CALC.(DOUBLE ASSAYS)**

**asc File Example:**

Definition: 96 rows  
 Column 1 Position  
 Column 2 sample name (or „empty“ if sample position not used)  
 Columns 3,4 not used  
 Columns 5-9 Data

A1	CAL 0.00		0.067	<Min	<Min		
B1	CAL 0.00		0.069	<Min			
C1	CAL 0.313		0.109	<Min	0.3331		6.4218
D1	CAL 0.313		0.121	0.3331			
E1	CAL 0.625		0.223	0.65434	0.62479	6.6885	-0.033407
F1	CAL 0.625		0.203	0.59524			
G1	CAL 1.25		0.449	1.271	1.2499	2.38	-0.0050019
H1	CAL 1.25		0.433	1.2289			
A2	CAL 2.50		0.958	2.474	2.5001	1.4778	0.0048336
B2	CAL 2.50		0.979	2.5262			
C2	CAL 5.00		1.663	5.0925	5.0005	2.6038	0.0094718
D2	CAL 5.00		1.625	4.9084			
E2	CAL 10.0		2.361	10.039	10	0.54452	0.0014085
F2	CAL 10.0		2.354	9.9616			
G2	CAL 20.0		2.759	>Max	18.827		-5.8664
H2	CAL 20.0		2.732	18.827			
A3	QCS 0.750		0.247	0.72357	0.71499	1.6976	-4.668
B3	QCS 0.750		0.241	0.70641			
C3	QCS 7.50		2.273	9.1401	8.6304	8.3535	15.071
D3	QCS 7.50		2.155	8.1206			
E3	QCS 15.0		2.683	16.26	16.599	2.8895	10.662
F3	QCS 15.0		2.699	16.938			
G3	empty						
H3	empty						
A4	QCS 0.750		0.234	0.68627	0.71771	6.1966	-4.3049
B4	QCS 0.750		0.256	0.74916			
C4	QCS 7.50		2.11	7.7737	8.3231	9.3361	10.975
D4	QCS 7.50		2.244	8.8726			
E4	QCS 15.0		2.688	16.461	16.942	4.0176	12.95
F4	QCS 15.0		2.709	17.424			
G4	empty						
H4	empty						
A5	QCS 0.750		0.255	0.74633	0.75622	1.8507	0.82971
B5	QCS 0.750		0.262	0.76612			
C5	QCS 7.50		2.113	7.7962	8.3798	9.8497	11.731
D5	QCS 7.50		2.254	8.9635			
E5	QCS 15.0		2.811	>Max	>Max		>Max
F5	QCS 15.0		2.821	>Max			

## **Appendix 7: Access2**

Access2 results are copied (pcAnywhere) from the Access2 system into pre-defined folder on the file server.

The \*.csv file is imported into dbLabCal via menu BATCH IMPORT...

***Important: csv file import is possible only if current project uses regression model  
NO CALC.(DOUBLE ASSAYS)***

## **Appendix 8: ISE**


The \*.xls file is imported into dbLabCal via menu BATCH IMPORT...

***Important: xls file import is possible only if current project uses regression model  
LOG.REGRESSION***

## Appendix 9: ICP-MS Elan

### Define Method (Icon „Method“ ).

Define the export rep file with „Report Filename“.



Define export settings in „Report Option Template“. The rop file saves the export settings.

Select following sections in the rop file „Report Option“:

- Sample Information
- Mean Values

Save report method e.g. as dlabelcalreport.rop (Report Template)

See following screen shots for details:

Report Options - C:\eladata\ReportOptions\EDV report.rop[Modified]

Available Sections		Selected Sections	
1	Sample Information	1	Sample Information
2	Replicates	2	Mean Values
3	Mean Values	3	
4	Standard Deviations	4	
5	Relative Std. Dev.	5	
6	Summary	6	
7	Calibration	7	
8	QC Calculated Values	8	
9	QC Out Of Limits	9	
10	QC Action	10	
11	Raw Data	11	
12	Dual Detector Gains		

Append -> Insert -> Remove

Available Fields		Selected Fields	
36	Number of Replicates	1	Batch ID
37	QC Mode	2	Sample ID
38	Dual Detector Mode	3	Sample Type
39	Blank subtracted	4	Sample Description
40	Measurement Unit	5	Acquisition Date/Time-Intl.
41	Peak Processing Mode	6	Sample File
42	Smooth Signal Profile Factor	7	Method File
43	Signal Profile Processing Mode	8	Aliquot Volume
44	Number of Baseline Readings	9	Diluted To Volume
45	Auto Lens Mode	10	Autosampler Position
46	User Name	11	User Name
47	Computer Name	12	Computer Name
48	Plain Text	13	
49	Report Date/Time	14	
50	Report Date/Time & Timezone	15	
51	Report Date/Time-Short	16	
52	Report Date/Time&Zone-Short	17	
53	Report Date/Time-Intl.	18	
54	Report Date/Time&Zone-Intl.	19	
55		20	
56		21	

Append -> Insert -> Remove

Report Options - C:\eladata\ReportOptions\EDV report.rop[Modified]

Available Sections		Selected Sections	
1	Sample Information	1	Sample Information
2	Replicates	2	Mean Values
3	Mean Values	3	
4	Standard Deviations	4	
5	Relative Std. Dev.	5	
6	Summary	6	
7	Calibration	7	
8	QC Calculated Values	8	
9	QC Out Of Limits	9	
10	QC Action	10	
11	Raw Data	11	
12	Dual Detector Gains		

Append -> Insert -> Remove

Available Fields		Selected Columns/Rows	
1	Internal Standard Symbol	1	Internal Standard Symbol
2	Analyte	2	Analyte
3	Mass	3	Mass
4	Measured Intensity Mean	4	Measured Intensity Mean
5	Blank Intensity	5	Blank Intensity
6	Net Intensity Mean	6	Net Intensity Mean
7	Concentration Mean	7	Concentration Mean
8	Report Unit	8	Report Unit
9	Mode	9	Mode
10		10	
11		11	
12		12	
13		13	
14		14	
15		15	
16		16	
17		17	
18		18	
19		19	
20		20	
21		21	

Append -> Insert -> Remove



### Write Sample list (Icon „Sample“).

Here are only the unknown samples. Standards and QCs are defined in the method.

ELAN Instrument Control Session - [Samples - C:\elandata\Sample\Lithium\SA061\SA061 seq 8.sam]

File Edit Analysis Options Wizard Window Help

Method Sample Dataset Realtime Interactive CalView RptOption RptView Optimize Tuning Instrument Devices SmartTune

Manual Batch

Analyze Batch Sample Template... Summary... Build Run List...

Batch Index	A/S Loc.	Batch ID	Sample ID	Measurement Action (*)	Method (*)	Description	Sample Type (*)
1	24		doubleblank	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
2	25		wash	Run Blank, Stds. and Sample	lithium\sa061\sa061 seq 8.mth		Sample
3	38		sub15dg01tp0h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
4	39		sub15dg01tp0.5h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
5	40		sub15dg01tp1h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
6	41		sub15dg01tp1.5h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
7	42		sub15dg01tp2h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
8	43		sub15dg01tp2.5h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
9	44		sub15dg01tp3h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
10	45		sub15dg01tp4h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
11	46		sub15dg01tp5h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
12	47		sub15dg01tp6h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
13	48		sub15dg01tp7h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
14	49		sub15dg01tp8h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
15	50		sub15dg01tp9h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
16	51		sub15dg01tp10h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
17	52		sub15dg01tp12h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
18	53		sub15dg01tp16h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
19	54		sub15dg01tp24h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
20	55		sub15dg01tp48h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
21	56		sub15dg01tp72h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
22	57		sub15dg01tp96h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample
23	58		sub15dg01tp120h	Run Sample	lithium\sa061\sa061 seq 8.mth		Sample

### Define the result file „Dataset“ (Icon „Dataset“)

ELAN Instrument Control Session - [Dataset - C:\elandata\DataSet\SA061 seq8]

File Edit Analysis Options Wizard Window Help


Method Sample Dataset Realtime Interactive CalView RptOption RptView Optimize Tuning Instrument Devices SmartTune

Reprocess Summary Report... Method Load


Use original conditions

Batch Index	Batch ID	Sample ID	Acquisition Date/Time	Method	Description	Read Type (*)	Sample File Name
1		doubleblank	13-May-2008 7:41:39 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	doubleblank.001
2		0 ng	13-May-2008 7:45:20 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Blank	0 ng.002
3		25 ng	13-May-2008 7:49:00 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #1	25 ng.003
4		50 ng	13-May-2008 7:52:41 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #2	50 ng.004
5		100 ng	13-May-2008 7:56:23 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #3	100 ng.005
6		500 ng	13-May-2008 8:00:05 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #4	500 ng.006
7		1000 ng	13-May-2008 8:03:45 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #5	1000 ng.007
8		2500 ng	13-May-2008 8:07:24 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #6	2500 ng.008
9		5000 ng	13-May-2008 8:11:03 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #7	5000 ng.009
10		10000 ng	13-May-2008 8:14:42 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Standard #8	10000 ng.010
11		7500 ng	13-May-2008 8:18:22 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		QC Std #1	7500 ng.011
12		wash	13-May-2008 8:22:03 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	wash.012
13		sub15dg01tp0h	13-May-2008 8:25:44 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp0h.013
14		sub15dg01tp0.5h	13-May-2008 8:29:24 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp0.5h.014
15		sub15dg01tp1h	13-May-2008 8:33:05 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp1h.015
16		sub15dg01tp1.5h	13-May-2008 8:36:45 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp1.5h.016
17		sub15dg01tp2h	13-May-2008 8:40:25 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp2h.017
18		sub15dg01tp2.5h	13-May-2008 8:44:05 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp2.5h.018
19		sub15dg01tp3h	13-May-2008 8:47:45 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp3h.019
20		sub15dg01tp4h	13-May-2008 8:51:25 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp4h.020
21		sub15dg01tp5h	13-May-2008 8:55:06 AM	C:\elandata\Method\lithium\sa061\sa061 seq 8.mth		Sample	sub15dg01tp5h.021

Select „Sample“ icon again, select all and click „Build Run List“ button



Print „Run List“



Click „Analyze Batch“button. **Analysis is running...**

All data are in the rep file after the analysis is finished.

The \*.rep file is imported into dbLabCal via menu BATCH IMPORT...

All data including the "MEAS. INTENS. MEAN" values as well as the "NET INTENS. MEAN" values are in a rep file after the measurement was finished.

dbLabCal imports either "MEAS. INTENS. MEAN" or "NET INTENS. MEAN" depending on the column position (the most left column data are imported) which can be defined while creating the rep file.

**Appendix 10: ICP-HPLC-MS Chromera**Designed for ONLY 1 Analyte plus 1 Internal Standard

In Excel „Search&Replace „ISName“ mit “ IS “+„ISName“

*dbLabCal identifies internal standard (one of the 2 allowed entries) if the element name starts with “IS “ characters vor*

Sample Name	Sample Description	Injection Number	Vial Number	Sample Type	Standard Name	Group
QCs 150 ng		1	1	Sample		UA29 0

Channel	Ret. Time	Component Name	Area	Height
Gd 160	17.14168111		129811.6075	1507.66624
IS Lu 175	19.32693663	Lutetium	1370701.828	10672.56676

QCs 1500 ng		1	2	Sample		UA29 0
-------------	--	---	---	--------	--	-----------

Channel	Ret. Time	Component Name	Area	Height
Gd 160	17.03649412		1236305.684	14069.79958
IS Lu 175	19.3152	Lutetium	1415110.746	10364.69024

QCs 4000 ng		1	3	Sample		UA29 0
-------------	--	---	---	--------	--	-----------

Channel	Ret. Time	Component Name	Area	Height
Gd 160	17.09205576		3432714.392	39252.48262
IS Lu 175	19.3431267	Lutetium	1515799.547	10645.02664

wash		1	4	Sample		UA29 0
------	--	---	---	--------	--	-----------

Channel	Ret. Time	Component Name	Area	Height
Gd 160	17.01013632		28490.2845	331.671417
IS Lu 175	19.17466486	Lutetium	6196.9763	61.08675937

wash		1	5	Sample		UA29 0
------	--	---	---	--------	--	-----------

Channel	Ret. Time	Component Name	Area	Height
Gd 160	17.12721705		12902.0261	157.6733076
IS Lu 175	19.45235138	Lutetium	3814.75719	33.17419784

wash		1	6	Sample		UA29 0
------	--	---	---	--------	--	-----------

USW...

## **Appendix 11: FACS**

Export results (program FlowCytomix Pro 2.3, view EVALUATION RESULTS) into a csv-file.

The \*.csv file is imported into dbLabCal via menu BATCH IMPORT...

### ***Important:***

***dbLabCal imports only calculated results. Therefore, xls file import is possible only if current project uses regression model No CALC.***

## **Appendix 12: SearchLight**

SearchLight software creates results in a xls file. xls file has also a special sheet adapted for dbLabCal import called PG-ML DATA2.

dbLabCal imports data from PG-ML DATA2 sheet automatically when selecting the xls file via dbLabCal's BATCH IMPORT... menu

dbLabCal checks also automatically double assay CV acceptance criterion if activated – see below- and set # flag if necessary.

	Standards	QC's
Max. deviation from nom. conc. at LOQ / LQC	20 %	15 %
Max. deviation from nom. conc. (remaining CALs/QC's)	15 %	15 %
Min. number of used and accepted values	6	4
		1
Set flag to # (all samples except CAL) if CV >	20 %	
Set flag for QC's/VALs to E if deviation higher than	off %	
Set flag for SUBs to S if IS not within mean(CAL) ±	off %	S
Set flag for QC's to S if IS not within mean(CAL) ±	off %	S

At least 1 QC's at the same concentration must fit acceptance criterion!

***Important: xls file import is possible only if current project uses regression model NO CALC.(DOUBLE ASSAYS)***

**Appendix 13: Mesoscale (MSD)**

Die MSD Discovery Bench Software creates one CSV file containing all results of all analytes. Analyte names are defined in the second column („Assay“). See following screen shot for an example of a MSD export file.

A2 Well										
A	B	C	D	E	F	G	H	I	J	
1	Plate *2CC2QAA038*									
Well	Assay	Sample Group	Sample	Concentration	Dilution	Signal	Mean	Std. Deviation	CV	
3	A03	TNF-? (Human)	Standard	S009	0	206	228	30.40559159	13.3650952	
4	A04	TNF-? (Human)	Standard	S009	0	249	228	30.40559159	13.3650952	
5	H01	TNF-? (Human)	Standard	S008	4	830	824	8.485281374	1.029767157	
6	H02	TNF-? (Human)	Standard	S008	4	818	824	8.485281374	1.029767157	
7	G01	TNF-? (Human)	Standard	S007	10	1418	1351	94.75230868	7.013494351	
8	G02	TNF-? (Human)	Standard	S007	10	1284	1351	94.75230868	7.013494351	

K	L	M	N	O	P	Q
Calc. Concentration	Calc. Conc. Mean	Calc. Conc. Std. Deviation	Calc. Conc. CV	% Recovery	% Recovery Mean	Detection Range
0	0.08504436	0.120270888	141.4213562			Below Fit Curve Range
0.170088721	0.08504436	0.120270888	141.4213562			Below Detection Range
5.391467711	5.342448938	0.069323014	1.297588706	134.7866928	133.5612234	In Detection Range
5.293430164	5.342448938	0.069323014	1.297588706	132.3357541	133.5612234	In Detection Range
10.03301317	9.514880521	0.732750223	7.701097465	100.3301317	95.14880521	In Detection Range
8.996747869	9.514880521	0.732750223	7.701097465	89.96747869	95.14880521	In Detection Range

dbLabCal creates individual batch for every imported analyte. This is the difference to the standard behaviour of dbLabCal: dbLabCal creates one batch containing all analytes analysed in one batch.

**Appendix 14: Calculation of weighted linear regression**

$$y = a + b \cdot x$$

---

 Weighting:
 

---

none:  $w_i = 1$

1/x:  $w_i = \frac{1}{x_i}$

1/y:  $w_i = \frac{1}{y_i}$

1/x<sup>2</sup>:  $w_i = \frac{1}{x_i^2}$

1/y<sup>2</sup>:  $w_i = \frac{1}{y_i^2}$ 


---

$$S_w = \sum_{i=1}^n w_i$$

$$S_x = \sum_{i=1}^n w_i \cdot x_i$$

$$S_y = \sum_{i=1}^n w_i \cdot y_i$$

$$S_{xy} = \sum_{i=1}^n w_i \cdot x_i \cdot y_i$$

$$S_{x^2} = \sum_{i=1}^n w_i \cdot x_i^2$$

$$S_{y^2} = \sum_{i=1}^n w_i \cdot y_i^2$$


---


Intercept: 
$$a = \frac{S_y \cdot S_{x^2} - S_x \cdot S_{xy}}{S_w \cdot S_{x^2} - S_x^2}$$

Slope: 
$$b = \frac{S_w \cdot S_{xy} - S_x \cdot S_y}{S_w \cdot S_{x^2} - S_x^2}$$

Determination Coefficient: 
$$r^2 = \frac{(S_w \cdot S_{xy} - S_x \cdot S_y)^2}{(S_w \cdot S_{x^2} - S_x^2) \cdot (S_w \cdot S_{y^2} - S_y^2)}$$


---

**Appendix 15: Flowchart: Reporting Final Results (Lang&Bolton)**





**Appendix 16: Flowchart: Reporting Final Results (HoLaRo)**