

Hardware and Software Requirements

dbLabCal V3 client runs on any computer which is suitable for current Windows operating system. Please refer also to the administration manual (dbLabCal_Administration.pdf)

	Minimum	Recommended
Processor	Intel (x86)	N/A
Memory	256 MB	2 GB
OS	Microsoft Windows XP Pro	W7 Pro or W8.1 Pro or higher, 64bit
Disk space	1.5 GB minimum	1.5 GB minimum
Windows Installer	MSI version 2.0 minimum	N/A
Network	TCP/IP	N/A
Resolution	1024x768	1280x1024 or higher
Other SW	MS Word and MS Excel are required for reporting and data import and should be installed on the client PC.	

Procedure

#	TASK	EXPECTED RESULT	ACTUAL RESULT
1	Record identification information for the client machine.	The data is recorded and all fields are complete.	PC Name: _____ IP Address: _____
2	Verify that the client machine fulfils minimum requirements	Minimum requirements fulfilled.	Processor: _____ Memory: _____ OS: _____ Disk space: _____ Resolution: _____ PASS FAIL
3	Install Oracle Client 11 (32bit!) recommended installation option: "Administrator" (from dbLabCal's setup folder)	Oracle client installed without errors/warnings.	 PASS FAIL
4	Install dbLabCal Client (from dbLabCal's setup folder)	dbLabCal client installed without errors/warnings.	 PASS FAIL
5	Create shortcut on desktop to dblabcal.exe in dbLabCal's bin folder on the server share	N/A	N/A
6	Start dbLabCal	dbLabCal starts with login dialog.	 PASS FAIL

Signature / Date