

Hardware and Software Requirements

dbLabCal V3 client runs on any computer which is suitable for current Windows operating system. Please refer also to the administration manual (dbLabCal_Administration.pdf)

	Minimum	Recommended
Processor	Intel (x86)	N/A
Memory	2 GB	16 GB
OS	Any W11 version or higher	Any W11 version or higher
Disk space	1.5 GB minimum	1.5 GB minimum
Resolution	1280x1024	1920x1200 or higher
Further software	MS Word (any version, but local installation!) is required for reporting functions. MS Excel might be required for data import from specific instruments	

Procedure

#	TASK	EXPECTED RESULT	ACTUAL RESULT
1	Record identification information for the client machine.	The data is recorded and all fields are complete.	PC Name: _____ IP Address: _____
2	Verify that the client machine fulfils minimum requirements	Minimum requirements fulfilled.	Processor: _____ Memory: _____ OS: _____ Disk space: _____ Resolution: _____ PASS FAIL
3	Install Oracle Client 11 (32bit!) recommended installation option: "Administrator" (from dbLabCal's setup folder)	Oracle client installed without errors/warnings.	PASS FAIL
4	Install dbLabCal Client (from dbLabCal's setup folder)	dbLabCal client installed without errors/warnings.	PASS FAIL
5	Create shortcut on desktop to dlabcal.exe in dbLabCal's bin folder on the server share <i>See: IQOQ dbLabCal Database, step 11</i>	N/A	N/A
6	Start dbLabCal	dbLabCal starts with login dialog.	PASS FAIL

Signature / Date